



## 9.7' LCD Pad Camera Microscope Accessories 12.0M A59.3520 With HDMI , Mini USB2.0 , USB 2.0 Interface

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.3520
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Theory: Digital Camera
- Resolution: 2048\*1536
- CPU: 1.8GHz Four Core RK3288 CPU
- WIFI: Support WIFI & Bluetooth
- Product Size: Pad 18.2\*23.7\*1.1cm
- Operation System: Android 5.1 System, Multi Language
- Highlight: **microscope accessory, compound microscope parts**



### More Images



## Product Description

### 9.7' LCD Pad Camera, 12.0M A59.3520 With HDMI, Mini USB2.0, USB 2.0 Interface


Specification	A59.3520 9.7' LCD Pad Camera
Digital Camera	12.0M, 2.2um*2.2um, Max 30 Fps
CPU	1.8GHz Four Core RK3288 CPU Support OpenGL ES 1.1/2.0 And OpenVG 1.1 Support Full HD 1080p Support Android 5.1 System
Memory	RAM 2GB, Support DDR3, DDR2, LPDDR2 ROM 8GB
Data Storage	SD Card Slot, Support Max 32GB C10 High Speed SD Card
Data Interface	HDMI, Mini USB2.0, USB 2.0
Operation System	Android 5.1 System, Multi Language
Software	<b>YX Camera Installed,</b> --Can Adjust Real Time Image --Can Measure Length, Area, Angle, Count --Can Take Picture, Record Video <b>YX Image Installed</b> --Grey, Contrast, Brightness Adjust --Binary Image For Granular Count  Both Software Chinese/English Version, With Instruction Manual Installed
Button	Power Switch
WIFI	Support WIFI, Support Bluetooth
Power	External Power Adaptor, Wide Voltage AC 100~240V, DC 5V/2A Power Mode : Normal, Sleep, Power Off
Product Size	Pad 18.2*23.7*1.1cm
Package Size	Packing Box 10*24*29cm, 1 pc/Box




A59.3520 9.7' LCD Pad Camera, 12.0M



A59.3520 9.7' LCD Pad Size



**A59.3520 9.7' LCD Pad Camera, 12.0M. Android System**  
 Built-in YW Camera Software, Can Take Photo, Capture Video, Measure Length, Angle And Etc.



**A59.3520 9.7' LCD Pad Camera, 12.0M. Android System**  
 Built-in YW Image Software, Can Adjust Grey, Contrast, Threshold to Measure Granular



A59.3520 9.7" LCD Pad Camera, 12.0M, Android System  
All Kinds Measuring Tools, Easy To Use



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China