



Student Optical Compound Microscope A11.1535 With LED Light Source , WF10X/18mm Eyepiece

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1535
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10X/18mm
- Certification: CE
- Drawtube: Monocular / Binocular / Trinocular
- Illumination: LED
- Theory: Biological Microscope
- Manification: 40x~1000x
- Highlight: **optical compound microscope,
wide field microscope**



Product Description

Student Biogocial Compound Optical Microscope A11.1535 With LED Light Source

A11.1535 is a very good choice for student, especially appreciated when budget is limited. All models are equipped with an integrated X-Y mechanical stage are supplied with an Abbe N.A. 1.25 condenser. Optional WF10X, WF16X, WF20X eyepiece enable utmost 2000X optical magnification, which can clearly observe biological slices, biological cells, bacteria and living tissue culture etc. And the curved handle is portable for users.



Parameter


A11.1535 Student Microscope				
Model	V	M	B	T
Head	Teaching Head	Monocular Head	Binocular Head	Trinocular Head
	Inclined at 30°, Interpupillary 48-75mm, Left Eye Tube Diopter Adjustable+/-5°			
Eyepiece	WF10X/18mm			
Nosepiece	Quadruple Nosepiece			
Objective	Achromatic Objectives 4x, 10x, 40x, 100x			
Condenser	Abbe Condenser NA1.25 with Diaphragm, with Filter Holder			
Focusing	Coaxial Coarse and Fine Adjustments, Moving Range 22mm			
Working Stage	Double Layer Mechanical Rectangular Stage 142X132mm, Moving Range 75X55mm			
Light Source	1W LED, Brightness Adjustable			
Adapter	1.0x C-Mount	N/A	N/A	1.0x C-Mount

Product Display

A11.1535-T Trinocular Head Student Microscope	A11.1535-M Monocular Head Student Microscope
A11.1535-B Binocular Head Student Microscope	A11.1535-V Teaching Head Student Microscope

SIZE AND WEIGHT:





PORTABLE:

 A11.1535-T Trinocular Head Student Microscope

Optional Accessories		
Eyepiece	WF16X/13mm	A51.1503-16
	WF20X/11mm	A51.1503-20
	WF10X/18mm Pointer Eyepiece	A51.1505
Objective	Achromatic Objectives 20x	A52.1505-20
	Achromatic Objectives 60x	A52.1505-60
	Plan Achromatic Objectives 4x, 10x, 40xs, 100xs(oil)	A52.1510
Light	6V/20W Halogm Lamp, Brightness Adjustable	A56.1504
	Upgrade to 3W LED, Brightness Adjustable	A56.1510
CCD Adapter	Eyepiece Tube Adapter	A55.1535-E
	0.5x	A55.1535-05

Plan Achromatic CF Objectives 4X 10X 40Xs 100Xs



A51.1503-16 A51.1503-20 A51.1505



A55.1535-E A55.1535-10 A55.1535-05



Product Details



01 HEAD

Monocular, binocular, trinocular and teaching head available for choose. All heads inclined at 30°, rotatable 360°, left eyepiece tube diopter adjustable +/-5°.

02 NOSEPIECE

Quadruple nosepiece with precise click stop. Equipped with 35mm achromatic CF objective 4X, 10X, 40Xs, 100Xs (oil).




03 WORKING STAGE

Double layer mechanical stage, size 142x132mm, moving range 75X55mm.

04 FOCUSING

Coaxial coarse & fine focusing adjustment, moving range 22mm.



05 LIGHT SOURCE

High brightness 1W LED, brightness adjustable.


06 POWER SWITCH & SOCKET

Available external power supply is AC100V-240V, frequency is 50/60Hz.



Product Application

Student microscope is designed for tertiary institutions and primary and secondary schools, which convenient to observe biological slices, biological cells, bacteria and living tissue culture, fluid precipitation observation and research. It is widely used in teaching demonstration, biochemical experiments, clinical research ect.



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China