



Monocular / Binocular Biological Microscope With Coaxial Coarse And Fine Focusing A11.1133

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1133
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student|School|Hobby|Home
- Certification: CE|Rohs
- Drawtube: Monocular/Binocular
- Illumination: LED
- Max Magnification: 1000X
- Theory: Biological Microscope
- Highlight: binocular biological microscope,
biological compound microscope



More Images



Product Description

Portable Monocular/Binocular Biological Microscope With Coaxial Coarse And Fine Focusing A11.1133

Specifications:

| A11.1133 Biological Microscope Specification | | A | B | C | D |
|--|---|---|---|---|---|
| Head | Monocular head, 30° inclined, 360° rotatable | ● | ● | | |
| | Sliding binocular head, 45° inclined, 360° rotatable | | | ● | |
| | Seidentopf binocular head, 30° inclined, 360° rotatable | | | | ● |
| Eyeiece | WF10X with pointer | ● | ● | ● | ● |
| Objective | Achromatic 4X | ● | ● | ● | ● |
| | Achromatic 10X | ● | ● | ● | ● |
| | Achromatic 40X, spring | ● | ● | ● | ● |
| | Achromatic 100X, spring/oil | | | ● | ● |
| Nosepiece | Triple revolving nosepiece | ● | | | |
| | Quadruple revolving nosepiece | | ● | ● | ● |
| Stage | Plain stage with clips, size 100x90mm | ● | | | |
| | Attachable mechanical stage, size 100x90mm | | ● | ● | |
| | Double layer Mechanical stage ,size 100x90mm | | | | ● |
| Condenser | Single lens condenser with Disc diaphragm | ● | | | |
| | N.A.1.25 condensor with Iris diaphragm | | ● | ● | ● |
| Focusing | Coaxial coarse and fine focusing knobs | ● | ● | ● | ● |
| Light Source | LED lampe with rechargeable battery | ● | ● | ● | ● |
| | Power supply adapter | ● | ● | ● | ● |



A11.1133-A





A11.1133-B

A11.1133-C



A11.1133-D



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China