



## View 1 - 95.2mm Stereo Binocular Microscopes A23.3645n-b1 Ce Rohs Certified

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.3645N-B1
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Eyepiece: WF10x/20mm
- Zoom Lens: 0.7~4.5x
- Zoom Ratio: 1:6.5
- Magnification: 7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives
- Working Distance: 100mm, With Auxiliary Objectives 30~165mm
- Field Of View: 1~95.2mm
- Stand: Pole Stand, No Light Source
- Highlight: **stereo zoom microscopes, stereo dissecting microscope**



### More Images




## Product Description

0.7-4.5x, The objective lens is doubled in ratio 6:4:1 to ensure image focus Binocular, WF10x/20, B1 Pole Stand + A1 Holder, Stand Size 200\*255\*22mm, Pole Height 240mm, Dia.32mm, No Light Source

| A23.3645N-B1 Zoom Stereo Microscope Specification |                                                                                                                               | B              | T |
|---------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------|----------------|---|
| Head                                              | Binocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm                         | O              |   |
|                                                   | Trinocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm, Eyepiece Tube Adapter |                | O |
| Eyepiece                                          | WF10x/20mm, Both With Eye-cup                                                                                                 |                |   |
| Zoom Lens                                         | 0.7~4.5x                                                                                                                      |                |   |
| Zoom Ratio                                        | 1:6.5                                                                                                                         |                |   |
| Magnification                                     | 7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Lens                                                            |                |   |
| Working Distance                                  | 100mm, With Standard 1x Lens                                                                                                  |                |   |
| Field Of View                                     | 1~95.2mm                                                                                                                      |                |   |
| Stand                                             | Pole Stand, No Light Source                                                                                                   |                |   |
| <b>Optional Accessories</b>                       |                                                                                                                               | Cata.No.       |   |
| Head                                              | Binocular Head + WF10x/20mm                                                                                                   | A53.3645N-B    |   |
|                                                   | Trinocular Head + WF10x/20mm                                                                                                  | A53.3645N-T    |   |
| Eyepiece                                          | WF10x/20mm, High-eyepoint                                                                                                     | A51.3621-1020  |   |
|                                                   | WF10x/22mm, High-eyepoint                                                                                                     | A51.3621-1022  |   |
|                                                   | WF10x/23mm, High-eyepoint                                                                                                     | A51.3621-1023  |   |
|                                                   | WF15x/15mm, High-eyepoint                                                                                                     | A51.3621-1515  |   |
|                                                   | WF20x/10mm, High-eyepoint                                                                                                     | A51.3621-2010  |   |
|                                                   | WF25x/9mm, High-eyepoint                                                                                                      | A51.3621-2509  |   |
|                                                   | WF30x/8mm, High-eyepoint                                                                                                      | A51.3621-3008  |   |
|                                                   | WF10x/22mm, High-eyepoint, Diopter Adjustable                                                                                 | A51.3622-1022  |   |
|                                                   | WF15x/16mm, High-eyepoint, Diopter Adjustable                                                                                 | A51.3622-1516  |   |
|                                                   | WF20x/12mm, High-eyepoint, Diopter Adjustable                                                                                 | A51.3622-2012  |   |
|                                                   | WF25x/9mm, High-eyepoint, Diopter Adjustable                                                                                  | A51.3622-2509  |   |
|                                                   | WF30x/8mm, High-eyepoint, Diopter Adjustable                                                                                  | A51.3622-3008  |   |
| Auxiliary Lens                                    | 0.5X, W.D. 165mm                                                                                                              | A52.3621-0.5x  |   |
|                                                   | 1.5X, W.D. 45mm                                                                                                               | A52.3621-1.5x  |   |
|                                                   | 2.0X, W.D. 30mm                                                                                                               | A52.3621-2.0x  |   |
| CCD Adapter                                       | 1.0x, For CTV                                                                                                                 | A55.3621-1.0x  |   |
|                                                   | 0.50x, For 1/2" CCD                                                                                                           | A55.3621-0.5x  |   |
|                                                   | 0.35x, For 1/3" CCD                                                                                                           | A55.3621-0.35x |   |
|                                                   | 0.50x, For 1/2" CCD                                                                                                           | A55.3621-N0.5x |   |

## PRODUCT DISPLAY



○ Eyepiece WF10X/20mm

○ 0.7x-4.5x, Magnification Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives



A54.3640-A1



A54.3630-B1

○ B1 Pole Stand + A1 Holder

**A23.3645N Zoom Stereo Microscope Optical Data Sheet**

| Objective | Standard 1.0x |            | Auxiliary 0.5X |            | Auxiliary 1.5X |            | Auxiliary 2X |            |
|-----------|---------------|------------|----------------|------------|----------------|------------|--------------|------------|
| W.D.      | 100mm         |            | 165mm          |            | 45mm           |            | 30mm         |            |
| Eyepiece  | Mag           | View Field | Mag            | View Field | Mag            | View Field | Mag          | View Field |
|           |               | 7.0X       | 28.6mm         | 3.5X       | 57.1mm         | 10.5X      | 19mm         | 14X        |

|          |       |        |        |        |         |        |      |        |
|----------|-------|--------|--------|--------|---------|--------|------|--------|
| 10X/20mm | 45.0X | 4.4mm  | 22.5X  | 8.9mm  | 67.5X   | 3mm    | 90X  | 2.2mm  |
| 15X/15mm | 10.5X | 21.1mm | 5.25X  | 42.8mm | 15.75X  | 14.3mm | 21X  | 10.7mm |
|          | 67.5X | 3.3mm  | 33.75X | 6.7mm  | 101.25X | 2.2mm  | 135X | 1.7mm  |
| 20X/10mm | 14.0X | 14.3mm | 7.0X   | 28.6mm | 21.0X   | 9.5mm  | 28X  | 7.1mm  |
|          | 90.0X | 2.2mm  | 45.0X  | 4.4mm  | 135.0X  | 1.5mm  | 180X | 1.1mm  |



**Opto-Edu (Beijing) Co., Ltd.**

0086 13911110627

[sale@optoedu.com](mailto:sale@optoedu.com)

[cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China