



## 45° Binocular Head Stereo Optical Microscope , 10x - 80x Stereo Binocular Microscope

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.3660N-B1
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Eyepiece: WF10x/20mm
- Head: 45° Binocular Head
- Objective: Turret Objective 1x3x Or 2X/4X
- Magnification Range: 10X-80X
- Working Distance: 100mm
- Stand: Pole Stand, No Light Source
- Highlight: **stereo dissecting microscope, stereo zoom microscopes**



A53.3660N Binocular Head + Body Magnification Up To 10x-80x.  
With WF10x/20mm Eyepiece

for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

Stereo microscope as a kind of visual equipment with three-dimensional, it plays an important role in biology, medicine, agriculture, forestry, industrial and marine biological industry

A22.3660N-B1 Stereo Microscope Specification	
Head	45° Binocular head
Eyepiece	Wide Field, High-eyepoint Eyepiece-WF10X20mm
Objective	Turret Objective 1x3x or 2X/4X
Magnification Range	10X-80X
Inter-Pupillary Distance	55-75mm
Diopter Range	5mm
Working Distance	100mm
Stand	Pole Stand, No Light Source

## PRODUCT DISPLAY



○ Eyepiece WF10X/20mm

- 20x/40x, Magnification Up To 10x.80x With Optional Eyepieces And 1X2X, 1X3X, 2X4X, Objectives



### Objective Optics Parameter

Model	Magnification	Eyepiece			Working Distance
		SZMEWh 10X20	SZMEWh 15X15	SZMEWh 20X10	
A22.3660N-B1-12	1X	10X20mm	15X15mm	20X10mm	100mm
	2X	10X20mm	30X7.5mm	40X5mm	
A22.3660N-B1-13	1X	10X20mm	15X15mm	20X10mm	
	3X	30X6.7mm	45X5mm	60X3.3mm	
A22.3660N-B1-24	2X	10X20mm	30X7.5mm	40X5mm	
	4X	40X5mm	60X3.75mm	80X2.5mm	



A54.3640-A1



A54.3630-B1



○ B1 Pole Stand + A1 Holder

A53.3660N-B



○ Binocular Head + Body Magnification Up To 10x.80x, With WF10x/20mm Eyepiece

A22.3660N-B1



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China