



## Eyepiece Stereo Optical Microscope Binocular Turret Objective Long Working Distance

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.3670-B1
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Eyepiece: WF10x/23mm
- Head: 45°Binocular Head
- Objective: Turret Objective 1x2x Or 1x3x Or 2X/4X
- Magnification Range: 10X-120X
- Working Distance: 110mm
- Stand: Pole Stand, No Light Source
- Highlight: **stereo dissecting microscope, stereo zoom microscopes**



## Product Description

The stereoscopic microscope has a positive stereoscopic sense, which is easy to operate and dissect. The field of view is large in diameter, the working distance is long, and the depth of focus is large, which is convenient to observe the whole layer of the object being detected.

A22.3670-B1 Stereo Microscope Specification	
Head	45° Binocular head
Eyeiece	Wide Field, High-eyepoint Eyepiece-WF10X23mm
Objective	Turret Objective 1x2x or 1x3x or 2X/4X
Magnification Range	10X-120X
Working Distance	110mm
Stand	Pole Stand, No Light Source



Objective Optics Parameter					
Magnification	Eyepiece				Working Distance
	10X23	10X20	15X16	20X12	
1X	10X23mm	10X20mm	15X16mm	20X12mm	110mm
2X	20X11.5mm	20X10mm	30X8mm	40X6mm	
3X	30X7.7mm	30X6.7mm	45X5.3mm	60X4mm	
4X	40X5.8mm	40X5mm	60X4mm	80X3mm	

## PRODUCT DISPLAY



○ Eyepiece WF10X/23mm

- 20x/40x, Magnification Up To 10x. 120x With Optional Eyepieces And 1X2X, 1X3X, 2X4X, Objectives





A54.3640-A1



A54.3630-B1



○ B1 Pole Stand + A1 Holder

A53.3670-B



○ Binocular Head + Body Magnification Up  
To 10x. 120x, With WF10x/23mm Eyepiece

A23.3670-B1



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China