



## 8 Stereo Microscope Accessories With 35 Stands Easy To Observe Objects

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A54.3630
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Use: 8 Stereo Microscope Series
- Quantity: 35 Stands
- Highlight: **compound microscope parts, microscope accessory**



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

8 Stereo Microscope Series + 35 Stands. The Head Can Be Combined Freely With Different Stands & Holder Arms. All Have Pole Stand, Pole Large Stand, Pole LED Stand, Pole Fan Stand, Track Stand, Track Large Stand, Track LED Stand, Track Fan Stand, Boom Stand, Single Arm Boom Stand, Dual Arm Boom Stand, Swing Arm Boom Stand.

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

**OPTO-EDU**

## 8 Stereo Microscope Series + 35 Stands

The head can be combined freely  
With Different Stands & Holder Arms



### Pole Stand



#### ○ **A54.3630-B1**

Stand 200\*255\*22mm, Pole Height 240mm,  
Dia.32mm, No Light Source.  
Used With Holder A54.3640-A1.

#### ○ **A54.3630-B2**

Stand 200\*255\*60mm, Pole Height 240mm,  
Dia.32mm, Up/Bottom 12V15W Halogen Lamp.  
Used With Holder A54.3640-A1.



#### ○ **A54.3630-B2L**

Stand 200\*255\*60mm, Pole Height 240mm,  
Dia.32mm, Up/Bottom 3W LED.  
Used With Holder A54.3640-A1.



#### ○ **A54.3630-B2LD**

Stand 200\*255\*60mm, Pole Height 240mm,  
Dia.32mm, Up 30 LED & Bottom 92 LED.  
Used With Holder A54.3640-A1.



#### ○ **A54.3630-B6**

Stand 200\*255\*22mm, Pole Height 240mm,



Stand 205x275x40mm, Pole Height 240mm,  
Dia.32mm, No Light Source.  
Used With Holder A54.3640-A1.

### ○ **A54.3630-B8LS**

Stand 205x275x40mm, Pole Height 240mm,  
Dia.32mm, Double Pipe LED.  
Used With Holder A54.3640-A1.



## Pole Large Stand



### ○ **A54.3630-B3**

Stand 380x260x20mm, Pole Height 240mm,  
Dia.32mm, No Light Source.  
Used With Holder A54.3640-A1.

### ○ **A54.3630-B3B**

Black Stand 380x260x20mm, Pole Height  
240mm, Dia.32mm, No Light Source.  
Used With Holder A54.3640-A1.



### ○ **A54.3630-B3BXLB**

Black Stand 380x260x20mm, Pole Height  
240mm, Dia.32mm, No Light Source.  
Used With Holder A54.3640-XLBA1.

### ○ **A54.3630-B3LS**

Stand 380x260x20mm, Pole Height 240mm,  
Dia.32mm, Double Pipe LED.  
Used With Holder A54.3640-A1.



## Pole LED Stand



### ○ **A54.3630-B8L**

Stand 205\*275\*40mm, Pole Height 240mm,  
Dia.32mm, Up/Bottom 3W LED.  
Used With Holder A54.3640-A1.

○ **A54.3630-B8LD**

Stand 205\*275\*40mm, Pole Height 240mm,  
Dia.32mm, Up 30 LED & Bottom 92 LED.  
Used With Holder A54.3640-A1.



○ **A54.3630-J4L**

Stand 285\*238\*40mm, Pole Height 240mm,  
Dia.32mm, Up/Bottom 3W LED.  
Used With Holder A54.3640-A1.

○ **A54.3630-J4LD**

Stand 285\*238\*40mm, Pole Height 240mm,  
Dia.32mm, Up 30 LED & Bottom 92 LED.  
Used With Holder A54.3640-A1.



Pole Fan Stand



○ **A54.3630-J1**

Stand 285\*238\*25mm, Pole Height 240mm,  
Dia.32mm, No Light Source.  
Used With Holder A54.3640-A1.

Track Stand

○ **A54.3630-B4**

Track Stand 200x255x22mm, Track Height  
300mm, Dia.76mm, Focusing Range 106mm,  
No Light Source.



○ **A54.3630-B7**

Track Stand 205x275x22mm, Track Height  
300mm, Dia.76mm, Focusing Range 106mm,  
No Light Source.

○ **A54.3630-B9LS**



Stand 205x275x40mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, Double Pipe LED.



○ **A54.3630-L**

Stand 255\*200\*28mm, Track Height 276mm, Dia.76mm, Focusing Range 78mm, No Light Source.



Track Large Stand

○ **A54.3630-B5**

Stand 380x260x20mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, No Light Source.



Track LED Stand



○ **A54.3630-B9L**

Stand 205\*275\*40mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, Up/Bottom 3W LED.

○ **A54.3630-B9LD**

Stand 205\*275\*40mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, Up 30 LED & Bottom 92 LED.



○ **A54.3630-J3L**

Stand 285\*238\*40mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, Up/Bottom 3W LED.

○ **A54.3630-J3LD**

Stand 285\*238\*40mm, Track Height 300mm, Dia.76mm, Focusing Range 106mm, Up 30 LED & Bottom 92 LED.





### ○ **A54.3630-2L**

Stand 255\*200\*40mm, Track Height 276mm,  
Dia.76mm, Focusing Range 74mm. LED Light.

## Track Fan Stand

### ○ **A54.3630-J2**

Stand 285\*238\*25mm, Track Height 300mm,  
Dia.76mm, Focusing Range 106mm, No Light  
Source.



## Boom Stand



### ○ **A54.3630-STL8**

Stand 380\*280\*10mm, Pole Height 355mm,  
Lateral movement 265mm, Dia.32mm, No Light  
Source.  
Used With Holder A54.3640-A3.

## Single Arm Boom Stand

### ○ **A54.3630-STL1**

Stand 230\*230\*40mm, Pole Height 380mm,  
Lateral movement 235mm, Dia.32mm, No Light  
Source.  
Used With Holder A54.3640-A1.



### ○ **A54.3630-STL3**

Stand 250\*240\*12mm, Pole Height 350mm,  
Lateral movement 400mm, Dia.38mm, No Light  
Source.  
Used With Holder A54.3640-A3.

### ○ **A54.3630-STL4**

Stand 260\*280\*15mm, Pole Height 380mm,  
Lateral movement 420mm, Dia.37mm, No Light  
Source.  
Used With Holder A54.3640-A3.



## Dual Arm Boom Stand



### ○ A54.3630-STL2

Stand 210\*260\*50mm, Pole Height 380mm,  
Lateral movement 235mm, Dia.32mm, No Light  
Source.  
Used With Holder A54.3640-A1.

### ○ A54.3630-STL5

Stand 260\*280\*15mm, Pole Height 380mm,  
Lateral movement 420mm, Dia.37mm, No Light  
Source.  
Used With Holder A54.3640-A3.



## Swing Arm Boom Stand



### ○ A54.3630-STL6A

Stand 98\*80\*116mm, vertical Height 400mm,  
Horizontal distance 770mm, Max clamping  
50mm, No Light Source.  
Used With Holder A54.3640-A3.

### ○ A54.3630-STL6B

Stand 98\*80\*116mm, vertical Height 400mm,  
Horizontal distance 770mm, Max clamping  
50mm, Dia.25mm, Pole Height 275mm, No  
Light Source.  
Used With Holder A54.3640-A3.



### ○ A54.3630-STL6C

Stand 260\*280\*15mm, vertical Height 400mm,  
Horizontal distance 770mm, Dia.25mm, Pole  
Height 323mm, No Light Source.  
Used With Holder A54.3640-A3.



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China