



7 - 45x Zoom Lens Stereo Optical Microscope With Up / Bottom 3w Led A23.3645-b2I

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.3645-B2L
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10x/20mm
- Zoom Lens: 0.7~4.5x
- Zoom Ratio: 1:6.5
- Magnification: 7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives
- Working Distance: 100mm, With Auxiliary Objectives 30~165mm
- Field Of View: 1~95.2mm
- Stand: Pole Stand, With Up/Bottom 3W LED
- Highlight: **stereo dissecting microscope, stereo zoom microscopes**



More Images



A54.3630-B2L Pole Stand + A54.3640-A1 Holder Arm

Product Description

0.7~4.5x, Binocular, WF10x/20, B2L Pole Stand + A1 Holder, Stand Size 200*255*60mm, Pole Height 240mm, Dia.32mm, Up/Bottom 3W LED, LED Source Low glare, no radiation, no pollution, no mercury, can be safely touched, a typical green lighting source, long life, can achieve colorful dynamic change effects and various images.



A23.3645-B2L Zoom Stereo Microscope Specification		B	T
Head	Binocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm	O	
	Trinocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm, Eyepiece Tube Adapter		O
Eyepiece	WF10x/20mm		
Zoom Lens	0.7~4.5x		
Zoom Ratio	1:6.5		
Magnification	7x-45x, Up To 3.5x-270x With Optional Eyepieces And Auxiliary Objectives		
Working Distance	100mm, With Auxiliary Objectives 30~165mm		
Field Of View	1~95.2mm		
Stand	Pole LED Stand, With Up/Bottom 3W LED		
Optional Accessories		Cata.No.	
Head	Binocular Head + WF10x/20mm	A53.3645-B	
	Trinocular Head + WF10x/20mm	A53.3645-T	
Eyepiece	WF10x/20mm, High-eyepoint	A51.3621-1020	
	WF10x/22mm, High-eyepoint	A51.3621-1022	
	WF10x/23mm, High-eyepoint	A51.3621-1023	
	WF15x/15mm, High-eyepoint	A51.3621-1515	
	WF20x/10mm, High-eyepoint	A51.3621-2010	
	WF25x/9mm, High-eyepoint	A51.3621-2509	
	WF30x/8mm, High-eyepoint	A51.3621-3008	
	WF10x/22mm, High-eyepoint, Diopter Adjustable	A51.3622-1022	
	WF15x/15mm, High-eyepoint, Diopter Adjustable	A51.3622-1515	
	WF20x/10mm, High-eyepoint, Diopter Adjustable	A51.3622-2010	
	WF25x/9mm, High-eyepoint, Diopter Adjustable	A51.3622-2509	
Auxiliary Lens	0.5X,W.D. 165mm	A52.3621-0.5x	
	1.5X,W.D. 45mm	A52.3621-1.5x	
	2.0X,W.D. 30mm	A52.3621-2.0x	
	1.0x, For CTV	A55.3621-1.0x	
	0.50x, For 1/2" CCD	A55.3621-0.5x	

CCD Adapter	0.35x, For 1/3" CCD	A55.3621-0.35x
	0.50x, For 1/2" CCD	A55.3621-N0.5x

PRODUCT DISPLAY



○ Eyepiece WF10X/20mm

- 0.7x-4.5x, Magnification Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives



A54.3640-A1



A54.3630-B2L



○ B2L Pole Stand + A1 Holder



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China