



High Precision Monocular Light Microscope Mechanical Stage 40X-1600X Magnification

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Model Number: A11.0106
- Minimum Order Quantity: 1
- Price: USD 672.5-433.9



Product Specification

- Magnification: 40X-1600X
- Head: Binocular Or Trinocular
- Eyepiece: WF 10X/18mm
- Nosepiece: Quadruple Revolving Nosepiece
- Stage: Mechanical Stage 150mmx135mm, Vernier Scale Travel Range 77mmx50mm With Graduation Of 0.1mm
- Highlight: **binocular compound microscope, wide field microscope**

for more products please visit us on cnoec.com

Product Description

1		
A11.0106	Specification	ABCD
Mag	nific 40X-1600X	• • • •
ation	Sliding Binocular Head, 45°Inclined, 360°Rotatable, Interpupillary Distance 55-75mm	• - - -
	Seidentopf Binocular Head, 45°Inclined, 360°Rotatable, Interpupillary Distance 55-75mm	- • - -
Hea	d Sliding Trinocular Head, 45°Inclined, 360°Rotatable, Interpupillary Distance 55-75mm	- - • -
	Seidentopf Trinocular Head, 45°Inclined, 360°Rotatable, Interpupillary Distance 55-75mm	- - - •
Eyep	iece WF10X/18mm	• • • •
	Achromatic 4X/0.10	• • - -
	Achromatic 10X/0.25	• • - -
	Achromatic 40X/0.65, Spring	• • - -
Obje	ctive Achromatic 100X/1.25, Spring, Oil	• • - -
	Semi-plan 4X/0.10	- - • •
	Semi-plan 10X/0.25	- - • •
	Semi-plan 40X/0.65, Spring	- - • •
	Semi-plan 100X/1.25, Spring, Oil	- - • •
Nos	epie Quadruple Revolving Nosepiece	• • • •
ce	Stag Mechanical Stage 150mm x 135mm, Vernier Scale Travel Range 77mm x 50mm with Graduation of 0.1mm,	• • • •
Con	dens Abbe N.A.=1.25, Milled Out to Phase Slide, Built-in Iris Diaphragm and Filter Holder, Vertical Rack and Pinion Focusing with Safely	• • • •
er	Spot, Height Adjustable	• • • •
	Blue, Amber and Green Filter	• • • •
Focu	sing Coaxial Coarse and Fine Focusing Knobs, 0.2mm per Rotation with Graduations in 0.002mm Increments, Tension Adjustment	• • • •
	Available, with Upper Stopper and Sliding Stopper	• • • •
	Light Base with Built-in Halogen Lamp 6V/20W , with Continue Intensity Sour Control, 110V/220V	• • • •
	ce	

Optional Accessories		Item No.	
Eyepiece	WF16X/10mm	EP16a	
	WF10X/18mm with Pointer	EPg10	
Objective	Reticle WF10X/18mm	EP10b	
	Plan 4X/0.10	B-OPr4	
	Plan 10X/0.25	B-OPr10	
	Plan 25X/0.4	B-OPr25	
	Plan 40X/0.65, Spring	B-OPr40	
Photo Attachment	Plan 100X/1.25, Spring, Oil	B-OPor100	
	Photo Adapter	MP03	
	Photo Eyepiece 2.5X	EPHn2.5	
CCD	Photo Eyepiece 4X	EPHn4	
	Photo Eyepiece 6.3X	EPHn6.3	
Dark Field	CCD C-Mount	TV02	
	Dark Field Condenser N.A.=0.83-0.91	C07	
Phase Contrast Attachment	Simple Dark Field Stop	DS03	
	Centering Telescope	EC8	
	Phase Contrast Condenser for Semi-Plan Phase Contrast Objective	C06a	
		10X/0.25	B-OSPnr10
	S-Plan Negative Phase Contrast Objective	40X/0.65	B-OSPnr40
		100x/1.25(Oil)	B-OSPnr100
		10X/0.25	B-OSPpr10
	S-Plan Positive Phase Contrast Objective	40X/0.65	B-OSPpr40
		100x/1.25(Oil)	B-OSPpr100
	Phase Contrast Condenser for Plan Phase Contrast Objective		C06
		10X/0.25	B-OPnr10
	Plan Negative Phase Contrast Objective	40X/0.65	B-OPnr40
		100x/1.25(Oil)	B-OPnr100
	Plan Positive Phase Contrast Objective	40X/0.65	B-OPpr10
		100x/1.25(Oil)	B-OPpr40
		OPopr100	



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China