



## OPTO-EDU A59.4249 Industrial HD Digital Camera HDMI Full HD Storage

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.4249
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Sensor: 1/2.86" Color CMOS
- Pixel Size: 2.75 $\mu$ m\*2.75 $\mu$ m
- Max Resolution: 1920x1080, 2.0M Full HD Preview, 8.0M High Resolution Photo Capture
- Definition: FULL HD
- White Balance: Manual Operation
- Memory Card: TF Class10
- Highlight: microscope accessory, parts of the microscope



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

A59.4249 HD Digital Camera Is HDMI Full HD Storage Industrial Camera, A High-speed 2.0 Million Pixel Preview And 8.0M Photo Resolution Color Digital Camera. And Large-format Array Cmos Image Sensor, Which Can Be Widely Used In Portrait And Human Biological Collection, Pattern Recognition, Fingerprint And Palm Print Collection, Industrial Detection, Medical Microscope Image And Machine Vision And Other Fields.

This Camera Include A Remote Controller To Take Photo, Record Video or Adjust Settings of Camera. , Which Helps To Avoid Shaking During Shooting Or Recording, The Images Captured Were Very Clear And Sharp .

The Camera Has Storage TF Slot, Can Use TF Class10 Or Above Memory Card, It Support Maximum 64GB TF Card. The Product Is Equipped With An Hdmi Adapter, Which Can Output Images Directly On Hdmi Monitor, TV, Projector, And Other Display Devices.Suitable For Teaching Operation Demon

Product Model	A59.4249
Sensor	1/2.86" Color CMOS
Pixel Size	2.75μm*2.75μm
Max Resolution	1920x1080, 2.0M Full HD Preview, 8.0M High Resolution Photo Capture
Definition	FULL HD
White Balance	Manual Operation
Memory Card	TF Class10
Photo/Video	Support
Cross Hair	Support
Menu	Brightness, Contrast, Saturation, RGB, Language Settings
Interface	HDMI/USB2.0
Power	DV-5V
Operating Temperature	-10°C+50°C
Operating Humidity	10%~90%

OPTO-EDU

# HD Digital Microscope Camera

A59.4249



- (●) Full HD 1080p Real Time Preview
- 📷 8.0M High Resolution Photo Capture
- ⚙️ TF Card Slot Support Photo & Video Save Easily
- 🔫 Infrared Controller To Take Photo & Adjust Parameter

Don't Miss Newest Camera For You!

# A59.4249 HD DIGITAL CAMERA

TF Card Slot


HDMI Image Output Interface

USB Power Cord Interface



Infrared Remote Control


1/2.86" Color CMOS




A59.4249 HDMI Camera + A23.1502 Zoom Stereo Microscope



A59.4249 HDMI Camera + A12.1502 Biological Microscope



A59.4249 HDMI Camera + A12.2601 Biological Microscope



Application Scenarios	
Biological Recognition	Portrait Detection, Fingerprint, Palmprint Recognition, Iris Detection, Etc
Defect Elimination	Test The Quality Of The Products
Accurate Measurement	Area Measurement, Three-coordinate Measurement, Three-dimensional Scanning, Stereo Vision Measurement, Etc
Tracking Identification	Qr Code, Commodity Code, Etc
Medical Research	Pathological Imaging, Microscopic Imaging, Specimen Collection, Etc

## PRODUCT DESCRIPTION



### How To Use



## PACKING AND SHIPPING

**Packing :** Cardbox And Polyfoam + Strong Master Carton 18\*13\*7cm, G.W.: 0.24Kgs  
**Shipping Time :** 3-10 Days From Stock



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China