




Binocular Optical School Student Microscope Monocular Trinocular Educational Gift A11.2601-B

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.2601
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Total Magnification: 40X 1000X
- Eyepiece: WF10X
- Head: Compensation Free Binocular, 30° Inclined, 360° Rotatable
- Objectives: 4x, 10x, 40x, 100x
- Nosepiece: Quadruple
- Focusing: Coaxial Coarse And Fine Focus Adjustment, Range: 22mm
- Stage: Double Layer Mechanical Stage 125mmx115mm, Moving Range : 75mmx35mm
- Condenser: ABBE NA1.25 Condenser
- Illuminator: 3W LED
- Highlight: **high powered microscope, binocular biological microscope, biological compound microscope**



More Images



for more products please visit us on cnoec.com

Product Description

Monocular Trinocular Educational Gift Cheap Binocular Optical School Student Microscope
Wide Field Eyepiece WF10x/18mm
High Contrast Achromatic Objective 4x,10x,40x,100x
Mechanical Stage 115x125mm, Moving Range 70x35mm
Full Metal Body High Rigid Construction
3W LED, Cordless Rechargeable Function Optional

A11.2601 Biological Microscope Is
Designed For Teaching, Clinical Verify.
The Series Microscope With Original Style,
Steady Structure, Convenient Operation
And Clear Image Is Suitable For Observing
Various Biological Specimens, And Mostly
Applied In School, College, Hospital And Family.

OPTO-EDU has our own factory for high end laboratory, polarizing, metallurgical, fluorescence microscope CNOPTec series, for comparison microscope A18 series, for SEM microscope A63 series, and for digital camera & LCD camera. Meantime, we are one of the best supplier for all kinds & levels microscope & accessories made in China, working in this field for more than 25 years, we have more than 1200 microscope models, and 5000 educational instrument to serve our customer with full range products.



WF10x-18mm, Eyepiece Is Anti-fungal
Anti-fungal System Built In Head

Head

WF10x-18mm, Eyepiece Is Anti-fungal, Anti-fungal System Built In Head



Achromatic Objective

Condenser N.A. 1.25

Mechanical Stage: 125mm × 115mm

Objective

Achromatic Objective: 4x, 10x, 40x(s), 100x/1.25(oil)(s), Objective Is Anti-fungal

Working Stage

Mechanical Stage, Stage Size: 125mm × 115mm, Travel: 35mm × 75mm

Abbe Condenser

N.A. 1.25, Adjustable Aperture, Aperture Center Is Adjustable, Built In Blue Filter, ϕ 38.5mm



Coaxial Coarse & Fine Focus

LED Illumination Systems

Focus

Coaxial Coarse & Fine Focus, Stage Moving Knob Can Move Stage Turn Left Or Right

Light Source

LED Illumination Systems (Non-rechargeable), Rotation Potentiometer With Power Switch



**A11.2601 Student Biological Microscope Upgrade To
A13.2603 Metallurgical Microscope**

A11.2601 Student Biological Microscope Upgrade To **A13.2603 Metallurgical Microscope**
A13.2603 Metallographic microscopes are often used to observe opaque objects such as wafers



**Powered By Safety Low Voltage Charger, Input Wide
Voltage 100V-240V, Output 5V1A, Support Power
Bank Supply For Outdoor Use**

Powered By Safety Low Voltage Charger

Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use

Items	A11.2601 Student Biological Microscope Specification	-M	-B	-T	-BP	-TP	Cata. No.
Head	Monocular Head Inclined 30°, Rotatable 360°	•					A53.2601-M
	Seidentopf Binocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm		•		•		A53.2601-B
	Seidentopf Trinocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm, Light Split Fixed P20/E80			•		•	A53.2601-T
Eyepiece	WF10x/18mm	•	••	••	••	••	A51.2602-1018
	WF16x/13mm	○	○	○	○	○	A51.2602-1613
	WF10x/18mm With Pionter	○	○	○	○	○	A51.2604-1018
	WF10x/18mm, With Micrometer X-0.1mm*180 Div / Y-Single Line	○	○	○	○	○	A51.2605-1018
	WF10x/18mm With Micrometer X-0.1mm*180 Div / Y-Single Line, Diopter Adjustable	○	○	○	○	○	A51.2605-1018D
Nosepiece	Quadruple, Backward	•	•	•	•	•	A54.2610-N03
Achromatic Objective	4x/0.1, W.D.=18.6mm	•	•	•			A52.2601-4
	10x/0.25, W.D.=6.5mm	•	•	•			A52.2601-10
	20x/0.40, W.D.=1.75mm	○	○	○			A52.2601-20
	40x/0.6(S), W.D.=0.47mm	•	•	•			A52.2601-40
	60x/0.80(S), W.D.=0.10mm	○	○	○			A52.2601-60
	100x/1.25(S,O), W.D.=0.07mm	•	•	•			A52.2601-100
Plan Objective	4x/0.1, W.D.=14.5mm				•	•	A52.2603-4
	10x/0.25, W.D.=5.65mm				•	•	A52.2603-10
	20x/0.40, W.D.=5.65mm				○	○	A52.2603-20
	40x/0.6(S), W.D.=0.85mm				•	•	A52.2603-40
	100x/1.25(S,O), W.D.=0.07mm				•	•	A52.2603-100
Stage	Double Layer Mechanical Size 125×115mm, Travel Range 35×75mm	•	•	•	•	•	A54.2601-S03
Condenser	Abbe N.A. 1.25(Iris Diaphragm)	•	•	•	•	•	A56.2614-09
Focusing	Coaxial Coarse And Fine Focusing Knobs, Travel Range: 22mm, Scale: 2µm	•	•	•	•	•	
C-mount	1.0x (Focus Adjustable)			○		○	A55.2601-B10
	0.5x (Focus Adjustable)			○		○	A55.2601-B05
Light Source	3W LED Illumination, Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use	•	•	•	•	•	A56.2613-3W
	3W LED Illumination (Rechargeable), Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use	○	○	○	○	○	A56.2651-3WR
	12V 20W Halogen Lamp, Wide Voltage: 100V-240V	○	○	○	○	○	A56.2650-20W
	Plano-Concave Mirror	○	○	○	○	○	A56.2601-B
Dark Field	Dark Field Ring, For 4x~40x Objective	○	○	○	○	○	A5D.2611-B
Polarizing	Analyzer + Polarizer Set	○	○	○	○	○	A5P.2601-B
Package	Carton Size 337*242*417mm, 1 Unit/Ctn						
Weight	Gross Weight 7.0kgs, Net Weigth 5.5kgs						

Note: "•"In Table Is Standard outfits, "○" Is Optional Accessories.

Product Models



A11.2601-M



A11.2601-B, BP



A11.2601-T, TP


Product Accessories



A5D.2611-B
Dark Field Ring



A56.2601-B
Plano-concave Mirror



A5P.2601-B
Analyzer / Polarizer



A55.2601-B05
0.5x C-Mount (Focus Adjustable)



A55.2601-B10
1x C-Mount (Focus Adjustable)



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China