



## Xsz-107bn Binocular Microscope Laboratory View Teaching 40X - 1600X Magnification

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1522-D
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Total Magnification: 40X~1600X
- Eyepiece: WF10X, WF16X
- Head: Compensation Free Binocular, 45° Inclined, 360° Rotatable (White Head)
- Objectives: 4x, 10x, 40x, 100x
- Nosepiece: Quadruple
- Focusing: Coaxial Coarse And Fine Focus Adjustment, Range: 30mm
- Stage: Double Layer Mechanical Stage 142mmx132mm, Moving Range : 75mmx45mm
- Condenser: ABBE NA1.25 Condenser With Diaphragm And Filter
- Illuminator: Halogen Lamp 6V/20W Adjustment
- Advantage: 40X-1600X Binocular Cheapest Classic Xsz-107bn Biological Microscope
- Highlight: **high powered microscope**,



### More Images



## Product Description

### OPTO-EDU A11.1522-D 40X-1600X Binocular Cheapest Classic xsz-107bn Biological Microscope

This product exports 2 million US dollars to Bangladesh!

A11.1522 Compound Biological Microscope Model, Can Be Widely Used In Medical And Sanitary Establishments, Laboratories, Institutes, Network Of Agricultural Research, Colleges And Universities For Clinical Microscope Inspections And Used As Instructional Equipment. They Also Can Be Used In Research Work Such As Biology, Bacteriology, Cytology, Histology And Pharmaceutical Chemistry.

OPTO-EDU is a unique, long-term and excellent supplier of microscopes and educational instruments from China. We focus on the international market of microscopes and education instruments. At present our customers come from more than 80 countries, such as the United States, Britain, Russia, Canada, Germany, Denmark, Poland, Sweden, the United Arab Emirates, Qatar, Saudi Arabia, Egypt, Argentina, Korea, Thailand, Brazil and so on. OPTO-EDU as our trademark is well known in the professional field, our products ranked No.1 in Alibaba.com. We have been working tirelessly for the sale of high quality optical products and educational products made in China to the world market.

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

**OPTO-EDU**


## Compound Biological Laboratory Microscope

**A11.1522**

- ☎ Cheapest Classic Model XSZ-107BN
- 🔍 WF10x Eyepiece, Achromatic 4x 10x 40x100x Objectives
- ⚙ Double Layer Mechanical Stage 142x132mm
- 💡 Halogen Lamp 6V/20W Adjustment

**Don't Miss Best Microscope For You!**

## A11.1522 COMPOUND BIOLOGICAL LABORATORY MICROSCOPE





A11.1522-A



A11.1522-B




A11.1522-C



A11.1522-D

#### A11.1522 Biological Microscope



# PRODUCT DISPLAY



## 01 HEAD

Compensation Free Binocular, 45° inclined, 360° Rotatable (White Head)

## 02 OBJECTIVE

Quadruple Nosepiece, Achromatic 4x, 10x, 40x, 100x Objectives Included



## 03 WORKING STAGE

Double Layer Mechanical Stage:  
132mmx142mm,  
Moving Range : 75mmx45mm

## 04 ILLUMINATOR

Halogen Lamp 6V/20W Adjustment, Abbe Condenser With Iris Diaphragm, With Filter Holder, And Moving Up/Down By Adjust Knob. Center Adjustable By 2 Screws On Condenser Holder



## 05 FOCUSING

Coaxial Coarse And Fine Focus Adjustment,  
Range: 30mm

## PRODUCT DESCRIPTION





A11.1522 Specification		A	B	C	D
Total Magnification	40X~1600X	•	•	•	•
Eyepiece	WF10X	•	•	•	•
	WF16X	•	•	•	•
Head	Monocular,30° Inclined,360° Rotatable	•	-	-	-
	Sliding Binocular,45° Inclined,360° Rotatable(Black Head)	-	•	-	-
	Sliding Trinocular,45° Inclined,360° Rotatable	-	-	•	-
	Compensation Free Binocular,45° Inclined,360° Rotatable(White Head)	-	-	-	•
Objectives	4x	•	•	•	•
	10x	•	•	•	•
	40x	•	•	•	•
	100x	•	•	•	•
Nosepiece	Quadruple	•	•	•	•
Focusing	Coaxial Coarse And Fine Focus Adjustment, Range:30mm	•	•	•	•
Stage	Double Layer Mechanical Stage 142mmx132mm, Moving Range : 75mmx45mm	•	•	•	•
Condenser	ABBE NA1.25 Condenser With Diaphragm And Filter	•	•	•	•
Illuminator	Halogen Lamp 6V/20W Adjustment	•	•	•	•



## OUR WORKSHOP



## OUR ADVANTAGES



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China