



Stereo Microscope Binocular 0.7-4.5x

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A23.1502-B1
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10x
- Zoom Lens: 0.7-4.5x
- Working Stage: Dia.100mm
- Focusing: Coarse
- Head: Binocular
- Light Source: N/A
- Highlight: stereo zoom microscopes, stereo binocular microscopes, stereo dissecting microscope



More Images



Product Description

OPTO-EDU A23.1502-B1 IT Testing Use Binocular Head Zoom Anatomical Microscope

A23.1502 Zoom Stereo Microscope Specification		B1	T1
Head	Binocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm	o	
	Trinocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm, With Eyepiece Adatper + 1x C-Mount		o

Our Product

Eyeiece	WF10x/20mm, High Eyepoint
Zoom Lens	0.7~4.5x
Zoom Ratio	1:6.5
Magnification	7x-45x, Up to 3.5x-180x With Optional Eyepieces and Auxiliary Objectives
Working Distance	100mm, With Auxiliary Objectives 30~177mm
Field of View	1~95.2mm
Light Source	No Light Source
Stand	Pole Stand ST1, Height 248mm, Base 200*255*22mm

OPTO-EDU


Zoom Stereo Microscope

A23.1502

- 0.7x~4.5x Zoom Lens, Zoom Ration 1:6.5
- WF10x/20mm High Eyepoint Eyepiece
- Classical Design High Quality Optical System
- More Than 35 Stands Optional

Don't Miss Best Microscope For You!

A23.1502 ZOOM STEREO MICROSCOPE





 A23.1502-B1 Binocular Zoom Stereo Microscope



 A23.1502-T1 Trinocular Zoom Stereo Microscope

Practical Application:

A23.1502 is a type of stereo microscope whose price is a little bit higher than others, because it is a professional stereo microscope in teaching demonstration, biological engineering, IT testing and some other fields. With zoom lens 0.7-4.5x, you can get the continuous zoom image you need. A23.1502 has high eyepoint with WF10x/20mm eyepiece, and the field of view is 1~95.2mm.

Zoom stereo microscope is widely used in teaching demonstration, bio-engineering, IT industry testing and other fields. A23.1502 has a high resolution and large field of view clarity, so it is also good at electronics industry and precision machinery industry parts assembly and testing, agricultural seed inspection and so on.

The coin under the stereo microscope:



STRUTURE & NOMENCLATURE



01 BIBOCULAR HEAD

Binocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm

02 TRINOCULAR HEAD

Trinocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm, With Eyepiece Adatper & 1x C-Mount



03 ZOOM

Zoom Lens 0.7~4.5x, Zoom Ratio 1:6.5

04 FOCUSING

Focusing Range 50mm, Distance to Center 155mm



05 WORKING STAGE

Pole Stand
Stand 200*255*22mm, Pole Height 240mm,
Dia.32mm



06 EYEPIECE ADAPTER & CCD ADAPTER

Eyeiece Adatper & 1x C-Mount

OPTIONAL ACCESSORIES



A51.1520-15
WF15x/15mm



A51.1520-20
WF20x/10mm



A55.1520-05
0.5x CCD Adapter



A52.1520-05
0.5X,W.D. 177mm



A52.1520-75
0.75X,W.D. 120mm



A52.1520-20
2.0X,W.D. 30mm

A23.1502 Zoom Stereo Microscope Optional Accessories		Cata. No.
Eyepiece	WF15x/15mm, High-eyepoint	A51.1520-15
	WF20x/10mm, High-eyepoint	A51.1520-20
Auxiliary Lens	0.5X,W.D. 177mm	A52.1520-05
	0.75X,W.D. 120mm	A52.1520-75
	2.0X,W.D. 30mm	A52.1520-20
CCD Adapter	0.5x	A55.1520-05

Company Information

OPTO-EDU (Beijing) Co., Ltd. manufactures many primary microscope and high-end laboratory, polarizing, metallurgical, fluorescence microscope, comparison microscope A18 series, for SEM microscope A63 series, and for digital camera & LCD camera with well-known brand OPTO-EDU and CNOPTEC. Meantime, we also supply full kinds microscope & accessories made in China, with experience of working in this field for more than 25 years, we have had more than 1200 models microscopes, and 5000 educational instruments in our range to serve our customer. Look forward to working with you!

OUR STRENGTH


- ◆ 25+ Years Professional Experience Know Microscopes Best
- ◆ 150+ Microscope & Accessory Manufacturers From China
- ◆ 200+ Hot Sale & Newest Update Microscope Model Every Month
- ◆ 720+ Customers From All Over World Keep Rising
- ◆ 1200+ Microscope Models To Create One-Stop Purchase Platform
- ◆ OPTO-EDU As Famous Brand For Microscope In China & USA
- ◆ Top Seller For Microscope on Alibaba Recommend To 99% User
- ◆ Most Professional Sale Team Supply Best Microscope Solution



Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China