



## A10.1007 Dark Field Microscope Sliding Trinocular Quadruple Achromatic

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A10.1007
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Magnification: 40x~1000x
- Eyepiece: Wide Field WF10x/18mm
- Nosepiece: Quadruple Nosepiece
- Filter: Blue, Yellow, Green
- Light Source: Built-in Illumination 6V20W Halogen Lamp
- Adapter: C-Mount 0.45x
- Highlight: **binocular compound microscope , wide field microscope**



## Product Description

### A10.1007 Dark Field Microscope Sliding Trinocular Quadruple Achromatic

A10.1007 Dark Field Microscope		A	B
Magnification	40x~1000x		
Head	Sliding Binocular Head, 45° Inclined, With C-Mount	○	
	Sliding Trinocular Head, 45° Inclined, With C-Mount		○
Eyepiece	Wide Field WF10x/18mm		
Nosepiece	Quadruple Nosepiece		
Objective	Achromatic 4x/0.10, 10x/0.25, 40x/0.65, 60X/0.85),100xs/1.25(Oil)		
	Plan 100X with Iris Diaphragm for Oil Dark Field Observation		
Stage	Double Layer Mechanical Stage, 140*140mm, Moving Range 50*75mm		
Focusing	Low Position Coaxial Coarse & Fine Focusing. Fine Focusing Scale 0.002mm, With Tension Adjustment, With Up Limit		
Condenser	NA1.25 Bright Field Condenser with Iris Diaphragm		
	Dark Field Condenser, N.A. 0.7~0.9, Dry		
	Dark Field Condenser, N.A. 1.25~1.35, Immersion		
Filter	Blue, Yellow, Green		
Light Source	Built-in Illumination 6V20W Halogen Lamp		
Adapter	C-Mount 0.45x		○
<b>A10.1007 Dark Field Microscope Optional Accessories</b>			
Adapter	Photo Adapter, For Nikon or Cannon SLR Camera, 1.0x	A55.1001	
	C-Mount 1.0x	A55.1002-10	
Eyepiece	P16x/11mm	A51.1002-16	




**A10.1007-A Dark Field Microscope**






Sliding Binocular Head, Diopter Adjustable



Quardpule Nosepiece, Achromatic Objectives



Built-in Illumination, Halogen Lamp 6V/ 20W



A11.1007-B



Plan Achromatic 100X with Iris Diaphragm for Oil Dark Field Observation



NA 0.7-0.9 For Dry Dark Field Condenser



NA 1.2-1.36 For Oil Dark Field Condenser



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China