



## Forensic Digital Bullet Comparison Microscope OPTO-EDU A18.1808 2.7x - 255x

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1808
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Magnification: 2.7x~255x
- Eyepiece: WF10X/22mm & WF20X/14mm
- Separation Line: Separation Line Adjusting System, Can Adjust Separation Line Easily
- Bullet Holder: Big Bullet Holder Set
- Adapter: SLR Camera Adapter & CCD Adapter
- CCD: Built-in Digital CCD Built-in, 3.0M, Video/USB Output
- Highlight: **bullet comparison microscope, binocular optical microscope**



### More Images



## Product Description

Forensic Digital Bullet Comparison Microscope OPTO-EDU A18.1808 2.7x - 255x

### Features:

Designed as big comparison microscope, **A18.1808XZB-12** comparison microscope, including computer control system and built-in digital CCD, it is suitable for forensic laboratory research, also good for tax, print, coinage, archaeology and ect.

- Total magnification 2.7x -255x
- Zoom objective 0.7x-4.9x
- 150W cold light for temperature sensitive objective view
- Tiltable big working stage, with transmit light hole on bottom
- Motorized focusing system can focus quickly
- Separation line adjusting system, can adjust separation line easily
- Built-in digital CCD built-in, 3.0M, video/USB output
- Electrical elevating working table optional

#### A18.1808 Comparison Microscope Specification

Magnification	2.7x~255x	
Eyepiece	WF10X/22mm	
	WF20X/14mm	
Objective	Zoom objective 0.7~4.9x, zoom ratio 1:7	
	Auxiliary lens 0.3x	
	Auxiliary lens 2x	
Separation Line	Separation line adjusting system, can adjust separation line easily	
Bullet Holder	Big bullet holder set	
Working Stage	Tiltable big working stage, with transmit light hole on bottom	
	Motorized focusing system	
Light Source	50W Air blast cooling light source provide high brightness	
	150W cold light source, soft optical fiber Dia.12mm with positioning ruler	
	Coaxial light lens	
	Polarizing lens	
Adapter	SLR Camera Adapter	
	CCD Adapter	
CCD	Built-in digital CCD built-in, 3.0M, video/USB output	
Optional Accessories		
Bullet Holder	Motorized bullet holder set	Item No.
Light Source	5W LED Light Source	A54.18WS-E
	5 Wave Band Light Source (red,yellow,blue,green,white)	A56.18LG-5W
	Hard Optical Fiber,pair, for quick adjustment of light spot	A56.18LG-H
	Ultraviolet light lens	A56.18LG-U
Adapter	Digital camera adapter for Canon A650	A55.18DC
Software	Professional forensic image analysis software <b>Crime Image</b>	A18.1815
Working Table	Electrical elevating working table	A54.18ST



Objective  
Auxiliary lens 0.3x, 2x

Eyepiece  
WF10x/22mm, WF20x/14mm

Ultraviolet light lens



Big bullet holder set

Red/Green Filter



5W LED Light Source



5 Wave Band Light



Hard Optical Fiber



Motorized bullet holder



Transmit light lens



Coaxial light lens



Polarizing lens

### Company Information

OPTO-EDU is a unique, long-term and excellent supplier of microscopes and educational instruments from China. We focus on the international market of microscopes and education instruments. At present our customers come from more than 80 countries, such as the United States, Britain, Russia, Canada, Germany, Denmark, Poland, Sweden, the United Arab Emirates, Qatar, Saudi Arabia, Egypt, Argentina, Korea, Thailand, Brazil and so on. OPTO-EDU as our trademark is well known in the professional field, our products ranked No.1 in Alibaba.com. We have been working tirelessly for the sale of high quality optical products and educational products made in China to the world market.

### OUR WORKSHOP





## OUR ADVANTAGES



As one of the most professional and most dynamic suppliers of optical and educational products in China, we have been concentrating on this field for more than 10 years. We have created the most professional product database of microscope & educational instrument, with 5000+ products and 500+ professional suppliers in system. We can provide our customers with the most comprehensive range of products in the related field, from the primary entry products to professional-level solutions, to meet the different needs of customers in the field of scientific research, teaching, industry, agriculture and various industries.



## OUR SERVICE



Feedback Within 12 Hours



After Sale Service Assured



3 Years Warranty Time

**contact information**

Welcome to visit our website, welcome your inquiry!

Manager	Mr. Xin Huang
Mobile phone	(+86)13911110627
Email	sale@cnoec.com
Website	www.cnoec.com.cn



**Opto-Edu (Beijing) Co., Ltd.**

 0086 13911110627     sale@optoedu.com     cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China