



A11.1513 1200X Children Biological Microscope Educational Student Microscope Set

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1513
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Color: Gray
- Eyepiece: WF10x
- Objective: 10x,40x,90x
- Working Stage: Plain Stage With Clips
- Focusing: Focusing
- Light Source: LED, Mirror
- Power Supply: Cordless Battery
- Nosepiece: 3 Holes
- Highlight: **high powered microscope, binocular biological microscope, biological compound microscope**



More Images



for more products please visit us on cnoec.com

Product Description

OPTO-EDU A11.1513 1200X Children's Biological Metal Microscope gift set; Kids educational student microscope set

A11.1513 Microscope is a beginner compound microscope. It has a monocular viewing head with LED and mirror illumination, a built-in color-filter wheel, a forward-facing rotating turret with 100x, 400x, and 900x magnification, and a white metal frame. The monocular viewing head has LED and mirror illumination and a built-in color-filter wheel for viewing a variety of specimen types. Coaxial coarse focus eases use for young users and has a rack-and-pinion focus mechanism for precise focusing. The plain stage has stage clips that secure the slide or specimen in place during viewing. The white metal frame is durable and stain-resistant. The microscope comes with a 52-piece accessory kit with brine shrimp eggs and hatchery, and a hard-sided plastic carrying case. The LED light is powered by two AA batteries (included).



COMPOND STUDENT MICROSCOPE

- Magnification 100x, 400x & 900x
- Plastic Suitcase Gift Box With 20 Standard Accessories
- Coaxial Coarse Focus
- Brightness LED And Mirror Light Source



Product Display



SIZE AND WEIGHT

	A11.1513	Plastic Box
Product Size	14x9x23CM	35x42x12CM
Net Weight	0.6KG	



 A11.1513 Compound Student Microscope



Accessories

Standard Accessories	
Packing Box	Plastic Suitcase Packing Box
Experimentation Appliance	Shrimp Hatchery
	Brine Shrimp Eggs
	Specimen Slicer
	Petri Dish
	Plastic Tweezers
	Plastic Scalpel
	Plastic Spatula
	Needle
	Plastic Stirring Rod
	Graduated Cylinder
	Eosin
	Gum Media (bottle)
	Sea Salt
	Prepared Slides*5
	Blank Slides*7
	Slide Covers*7
	Statistical Slide Covers*7
Slide Labels*7	
Collecting Vials*2	
Spare LED Bulb	



Product Details



01 HEAD

Monocular head with detachable WF10X eyepiece.

02 NOSEPIECE

Streamline outward triple, equipped with objective 10x, 40x, 90x.



03 WORKING STAGE

Plain stage with clips, size 65x72mm.

04 CONDENSER

Single lens condenser, 8 holes disc diaphragm make it easy to adjust the light source.



05 POWER SUPPLY

2 pcs AA batteries(included), the battery is in the base.



06 ILLUMINATION


LED and mirror light.

Related Products



Product Application

Student microscope is designed for tertiary institutions and primary and secondary schools, which convenient to observe biological slices, biological cells, bacteria and living tissue culture, fluid precipitation observation and research. It is widely used in teaching demonstration, biochemical experiments, clinical research etc.



Packaging & Shipping

PACKING: Polyfoam protection and strong carton.

Package	5 Sets/Carton, 55x40x36CM
Gross Weight	11KG
Packing	Polyfoam Protection And Strong Carton
Shipping Time	5~20 Days After Receive The Deposit



SHIPPING TIME: 5-20 days after receive the deposit.



Company Information



📍 OPTO-EDU is you're looking for, isn't it?



<p>25+ Yeas Professional Experience Know Microscopes Best</p> <p>150+</p> <p>150+ Microscope & Accessory Manufacturers From</p>	<p>200+ Hot Sale & Newest Update Microscope Model Every Month</p> <p>720+</p> <p>720+ Customers From All Over World Keep Rising</p>	<p>1200+ Microscope Models To Creat One-Stop Purchase Platform</p> <p>1200+</p>	<p>OPTO-EDU As Famous Brand For Microscope In China & USA</p> <p>Top Seller For Microscope on Alibaba Recommend</p>	<p>Most Professional Sale Team Supply Best Microscope</p>
--	--	--	---	---

Our Team And Our Factory.

OPTO-EDU, as one of the most professional manufacturer and supplier of microscope in China, our sub-brand CNOPTEC series high end biological, laboratory, polarizing, metallurgical, fluorescence microscopes, CNCOMPARISON series forensic

microscope, A63 series SEM microscope, and .49 series digital camera, LCD camera are very popular in the world market.



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China