



40x - 1000x A12.1010 Binocular Biological Microscope Laboratory Biological Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1010
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Optical System: Infinity Optical System
- Magnification: 40x~1000x
- Eyepiece: Wide Field WF10x/22mm, Eyepiece Tube Dia.30mm
- Nosepiece: Backward Quadruple Nosepiece
- Objective: Infinity Plan Achromatic 4x/0.10, 10x/0.25, 40xs/0.65, 100xs/1.25(Oil)
- Stage: Double Layer Mechanical Stage, 216*150mm, Moving Range 54*78mm, Double Slide Holder, Low Position X/Y Coaxial Control Knob, With Travel Guide To Avoid Scratch
- Condenser: N.A.1.25 Abbe Condenser, Iris Diaphragm, With Filter Holder, Up/Down Range 25mm, Center Adjustable, With Filter Holder
- Filter: Green
- Light Source: Halogen Lamp 12V20W , Brightness Adjustable



More Images



Product Description

A12.1010 Laboratory Biological Microscope, Newly designed with comfortable operating with ergonomic structure, with better illumination, convenient replacement of attachment with sliding-in center adjustable condenser.

- Infinity Optical System
- Seidentopf Binocular Head, 30° Inclined, 360° Rotatable
- Infinity Achromatic Objective 4x, 10x, 40x, 100x
- Abbe Condenser N.A.1.25, Center Adjustable
- Halogen Lamp 12V20W, Brightness Adjustable
- Low Position Coaxial Coarse & Fine Focus System



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085


A12.1010 Laboratory Biological Microscope



A12.1010 Laboratory Biological Microscope



A12.1010-A
Binocular, Quadruple
Laboratory Biological Microscope




A12.1010-B
Trinocular, Quadruple
Laboratory Biological Microscope



Seidentopf Trinocular Head, 30° Inclined, 360° Rotatable, Interpupillary 48~75mm, Wide Field WF10x/20mm, Dia.30mm



Backward Quadruple Nosepiece, With Infinity Plan Achromatic 4x/0.10, 10x/0.25, 40x/0.65(S), 100x/1.25(S,Oil)



Double Layer Mechanical Stage, 216x150mm, Moving Range 54x78mm, Double Slide Holder, Low Position X/Y Coaxial Control Knob, With Travel Guide To Avoid Scratch



Halogen Lamp 12V/20W, Brightness Adjustable
Light Source Cover Removable For Easy Repair



You Can See The Change Of Objective Magnification Even From The Back



Focusing
Low Position Coaxial Coarse & Fine Focusing



Knob
Switch & Light Source Adjustment knob

Items	A12.1010 Laboratory Biological Microscope	B	T
Head	Seidentopf Binocular Head, 30° Inclined, 360°Rotatable, Interpupillary Distance 48-75mm	•	
	Seidentopf Trinocular Head, 30° Inclined, 360°Rotatable, Interpupillary Distance 48-75mm		•
Eyepiece	Wide Field WF10x/20mm, Eyepiece Tube Dia.30mm, 2 pcs	•	•
Nosepiece	Backward Quadruple Nosepiece	•	•
Optical System	Infinity Optical System	•	•
Objective	Infinity Plan Achromatic 4x/0.10, 10x/0.25, 40x/0.65(S), 100x/1.25(S,Oil)	•	•
Magification	40x-1000x	•	•
Stage	Double Layer Mechanical Stage, 216x150mm, Moving range 54x78mm, Double Slide Holder, Low Position X/Y Coaxial Control Knob, With Travel Guide To Avoid Scratch	•	•
Focusing	Low Position Coaxial Coarse Fine Focusing. Fine Focusing Scale 0.002mm, Moving Range 25mm, With Safety Stop	•	•
Condenser	N.A.1.25 Abbe Condense, Iris Diaphragm, With Filter Holder, Up/Down Range 25mm, Center Adjustable, With Filter Holder	•	•
Filter	Green	•	•
Light Source	Halogen Lamp 12V/20W, Brightness Adjustable	•	•
Note: "•" In Table Is Standard Outfits			
Items	A12.1010 Laboratory Biological Microscope Optional Accessories		
Eyepiece	Wide Field WF15x/16mm, Eyepiece Tube Dia.30mm	A51.1002-15B	
	Wide Field WF20x/12mm, Eyepiece Tube Dia.30mm	A51.1002-20B	


Objective	Infinity Plan Achromatic 20x	A52.1008-20
	Infinity Plan Achromatic 60x	A52.1008-60
Nosepiece	Backward Quintuple Nosepiece	A54.1001-5
Illumination	Kohler Illumination	A56.1005-5
Adapter	Photo Adapter	A55.1008
	C-Mount 0.4x	A55.1009-04
	C-Mount 0.5x	A55.1009-05
Dark Field	Dark Field Condense, Dry	A5D.1012-D
	Dark Field Condense, Immersion	A5D.1012-I
Polarizing	Simple Polarizing Kit	A5P.1005
Phase Contrast	Phase Contrast Set	A5C.1006
Fluorescent	EPI-Fluorescent Kit	A5F.1005

Company Information



📍 OPTO-EDU is you're looking for, isn't it?

OPTO-EDU cnoptec



🏢 Our Team And Our Factory.

OPTO-EDU, as one of the most professional manufacturer and supplier of microscope in China, our sub-brand CNOPTEC series high end biological, laboratory, polarizing, metallurgical, fluorescence microscopes, CNCOMPARISON series forensic microscope, A63 series SEM microscope, and .49 series digital camera, LCD camera are very popular in the world market.



😊 I'm OPTO-EDU, nice to meet you!



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China