



USB 3.0 Port Microscope Accessories A59.4212 Digital Camera 1.85*1.85 Um Pixel Size

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.4212
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month




Product Specification

- Sensor: 1/1.7 Sony CMOS
- Pixel Size: 1.85*1.85 Um
- Resolution: 4000*3000
- Shutter Mode: Rolling Shutter
- Data Bits: 12 Bit
- Output Color: Color
- Output Frame Rate: Bin(2x2):30FPS/ALL-pixel:15FPS
- Interface: USB 3.0 Port
- Highlight: microscope accessory, parts of the microscope, compound microscope parts



for more products please visit us on cnoec.com

More Images



Product Description


30FPS 1/1.7 Sony CMOS USB 3.0 Port A59.4212 USB3.0 Digital Camera

A59.4212 Digital Camera, USB3.0	
Sensor	1/1.7 Sony CMOS
Pixel Size	1.85*1.85 um
Resolution	4000*3000
Shutter Mode	Rolling Shutter
Data Bits	12 bit
Output Color	Color
Output Frame Rate	Bin(2x2):30FPS/ALL-pixel:15FPS
Interface	USB 3.0 Port

Product Description



A59.4212 Digital Camera, USB 3.0



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China