



A34.4160 Mini Handheld USB Digital Optical Microscope 25x - 600x Rohs Certification

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.4160
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Color Of Body: Titanium Gray
- Magnification: 25x-600X (Manual)
- Focus Range: 0mm ~ 40mm
- Still Image Resolution: 1600*1200, 1280*960, 640*480, 480*320
- Lens: Microscope Lens , M12
- Resolution: 2.0M (Default 0.3M)
- Frame Rate: 30 F/s Under 600 Lux Brightness
- Image Sensor: Color CMOS Sensor
- Controller: High Speed DSP(Driver Free Available)
- Highlight: **digital usb microscope,**
portable digital microscope



More Images




Product Description

A34.4160 Mini handheld USB Digital Optical Microscope 25x - 600x

A 34.4160 Mini handheld USB Digital Microscope	
Color of Body	Titanium gray
Magnification	25x-600X (Manual)
Focus Range	0mm ~ 40mm
Still Image Resolution	1600*1200, 1280*960, 640*480, 480*320
Lens	Microscope Lens , M12
Resolution	2.0M (Default 0.3M)
Frame Rate	30 f/s Under 600 Lux Brightness
Digital Zoom	5X Sequence Mode
Image Sensor	Color CMOS Sensor
Controller	High Speed DSP(Driver Free available)
Power Supply	USB Port (5V DC)
Interface	USB 2.0
Illumination	White Light 8pcs LED
Snap Shot	Software and Hardware
Measure Tool	Micro-Measure(It can measure Length and diameters, angle, the perimeter, area, etc) (Language:English,Chinese Simplified)
Driver Language	English,Chinese
Operation System	Windows XP , Vista, Windows 7, 32 bit and 64 bit
System Requirement	Pentium Computer with 700M Hz & Above , 20M HD Space CD ROM Driver, 64MB RAM , Direct X VGA Card
Accessories	1.5M USB cable, CD disk, Metal Fix Stand, color box, calibration ruler, user manual



A34.4160



Calibration Ruler



8pcs LED



25x-600x



Circuit Board



Money



Texture




Skin




Printing

Company Information



📍 OPTO-EDU is you're looking for, isn't it?

OPTO-EDU *cnoptec*



Our Team And Our Factory.

OPTO-EDU, as one of the most professional manufacturer and supplier of microscope in China, our sub-brand CNOPTec series high end biological, laboratory, polarizing, metallurgical, fluorescence microscopes, CNCOMPARISON series forensic microscope, A63 series SEM microscope, and .49 series digital camera, LCD camera are very popular in the world market.



😊 I'm OPTO-EDU, nice to meet you!



👍 Trust us, we will give you a guarantee.

Direct factory price!

3 years quality warranty!

Payment by T/T, PayPal, West Union to LC!

Ranked No.1 supplier for microscope on alibaba.com!

Special price support OPTO-EDU/CNOPTEC distributor!

Opto-Edu (Beijing) Co., Ltd.

Transaction Level

41 Transaction(s)(6months) Response Rate

1,400,000+ 100.0%

The Supplier Supports Trade Assurance-A Free Service That Protects Your Orders From Payment To Delivery.

Supplier's Trade Assurance Limit: US \$379,000



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China