



1W LED Binocular Stereoscopic Microscope A22.1309 1W 110 - 240V Power Supply

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1309
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10x
- Objective: 2x/4x
- Working Stage: Dia.95mm
- Focusing: Coarse
- Light Source: 1W LED 18LEDs
- Power Supply: 110~240V
- Highlight: stereo zoom microscopes,
stereo binocular microscopes,
stereo dissecting microscope



More Images



Product Description

A22.1309 LED 40x Binocular Stereo Microscope

A22.1309 is the basic style stereo microscope with newly designed built-in Up Coaxial LED Light. Stereo microscope is also known as dissecting microscope, it is widely used in biology, medicine, agriculture, forestry, industry and marine life etc. A22.1309 has top and bottom light source, up coaxial LEDs inside the Objective, and the bottom light can adjust brightness. Its body is not huge, the track stand is convenient for users.

Parameter

A22.1309 Stereo Microscope	
Head	Binocular Head, Inclined 45°, Interpupillary Distance 55-75mm, With A Diopter Adaptor On The Right Eyepiece Tube
Eyepiece	WF10x/18mm
Objective	2x/4x
Working Distance	100mm
Stage	Ø95mm, 1 Pc Frosted Glass Stage(For Bottom Illumination)
Focus	Coarse Focusing, Focusing Range 100mm
Stand	Track Stand, Base Size 25x21x5cm, Height 31cm
Illumination	Up Coaxial LED, Bottom 1W LED (Brightness Adjustable), Up Oblique LED Optional
Power Supply	110~240V Wide Range Power Supply



STEREO MICROSCOPE

- Magnification 20x~40x Available
- ☀ Up Coaxial LED And Bottom LED Brightness Adjustable
- 📦 Portable With Track Stand



Product Display

SIZE AND WEIGHT

Product Size	25x21x41CM
Net Weight	4.1KG



 A22.1309 Stereo Microscope



01 BINOCULAR HEAD

Inclined 45°, interpupillary distance 55-75mm, with a diopter adaptor on the right eyepiece tube & WF10x/18mm eyepiece.

02 OBJECTIVE

2x/4x, with coaxial LED inside.



03 TRACK STAND

Base size 25x21x5cm, height 31cm.

04 POWER SUPPLY

110~240V wide range power supply.



05 LIGHT SOURCE

Two Switch Buttons For Up Coaxial LED And Bottom LED.

Bottom LED Light, brightness adjustable, 1 piece frosted glass stage for bottom illumination.

Stage Ø95mm.




Related Products



Product Application

Stereo microscope, also known as "solid microscope" or "operation and dissecting microscope", is a kind of stereoscopic microscope. Stereo microscope is widely used in materials such as macroscopic surface observation, failure analysis, fracture analysis and other industries field. As a kind of visual equipment with three-dimensional, it plays an important role in biology, medicine, agriculture, forestry, industrial and marine biological industry.



Packaging & Shipping

PACKING: Polyfoam protection and strong carton.

Package	1 Sets/Carton, 34X25X41CM
Gross Weight	5KG
Packing	Polyfoam Protection And Strong Carton
Shipping Time	5~20 Days After Receive The Deposit



SHIPPING TIME: 5~20 days after receive the deposit.



Company Information

OPTO-EDU is a unique, long-term and excellent supplier of microscopes and educational instruments from China. We focus on the international market of microscopes and education instruments. At present our customers come from more than 80 countries, such as the United States, Britain, Russia, Canada, Germany, Denmark, Poland, Sweden, the United Arab Emirates, Qatar, Saudi Arabia, Egypt, Argentina, Korea, Thailand, Brazil and so on. OPTO-EDU as our trademark is well known in the professional field, our products ranked No.1 in Alibaba.com. We have been working tirelessly for the sale of high quality optical products and educational products made in China to the world market.

OUR WORKSHOP



OUR ADVANTAGES

As one of the most professional and most dynamic suppliers of optical and educational products in China, we have been concentrating on this field for more than 10 years. We have created the most professional product database of microscope & educational instrument, with 5000+ products and 500+ professional suppliers in system. We can provide our customers with the most comprehensive range of products in the related field, from the primary entry products to professional-level solutions, to meet the different needs of customers in the field of scientific research, teaching, industry, agriculture and various industries.



OUR SERVICE



Feedback Within 12 Hours



After Sale Service Assured



3 Years Warranty Time



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China