



0.2x - 4.9x 3.0M Compound Optical Microscope OPTO-EDU A18.1814 Multi Function Comparator

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1814
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- OPTO-EDU A18.1814 Multi 0.2x~4.9x Function Comparator:
- Working Distance: 75~90mm
- View Field: Dia.1.5~60mm
- CCD: Built-in Digital CCD Built-in, 3.0M, Video/USB Output
- Software: Professional Forensic Image Analysis Software Crime Image
- PC: 14.1" TFT Notebook, With Software Installed
- Highlight: **binocular compound microscope, wide field microscope**



More Images



Product Description

0.2x~4.9x 3.0M With Software OPTO-EDU A18.1814 Multi Function Comparator

| | |
|--|--|
| A18.1814 GSB-4 Multi Function Comparator Specification | |
| Magnification | 0.2x~4.9x |
| Objective | Zoom objective 0.7x~4.9x, zoom ratio 1:7 Auxiliary lens 0.3x |
| Working Distance | 75~90mm |
| View Field | Dia.1.5~60mm |
| Bullet Holder | Small bullet holder set Motorized bullet holder set |
| Working Stage | Titable working stage, for bullet and small object Large motorized working stage, for document and large object Flourscence transmit light stage |
| Light Source | Excite & receive light source wave length: 254nm (Ultraviolet), 385nm (Ultraviolet), 540nm (Yellow), 680nm (Red), 800nm (Infrared), 960nm (Infrared), 1100nm (Infrared) Cold light source, for temperature sensitive objects viewing Holegen light source Flourscence ring light source Flourscence transmit light stage |
| CCD | Built-in digital CCD built-in, 3.0M, video/USB output |
| Software | Professional forensic image analysis software Crime Image |
| PC | 14.1" TFT notebook, with software installed |
| A18.1814 GSB-4 Multi Function Comparator Packing List | |
| Main Body | 1 set |
| Auxiliary lens 0.3x | 1 pc |
| Magnification correcting lens | 1 pc |
| Small bullet holder set | 1 set |
| Motorized bullet holder set | 1 set |
| Titable working stage, for bullet and small object | 1 set |
| Large motorized working stage, for document and large object | 1 set |
| Flourscence transmit light stage | 1 set |
| Flourscence ring light source | 1 set |
| Video cable | 1 pc |
| Power cord | 1 pc |
| Holgen lamp | 2 pcs |
| 2A Fuse | 2 pcs |
| Lens cleaning paper & Brush | 1 set |
| Instruction Manual | 1 book |

Our Product Introduction

for more products please visit us on cnoec.com





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China