



Opto Edu A43.1903 Portable Colposcope 6x Binocular 21v 150w Halogen Lamp

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A43.1903
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 6X
- Interpupillary Distance: 50mm-80mm
- Diopter: $\pm 5D$
- Magnification Changer: 3-step Magnification Change: 0.6X, 1X, 1.6X
- Focal Length Of Objective: F=200 & F=300mm (M45x0.75mm)
- Total Magnification: 3X, 5X, 8X, 4.7X, 7.5X, 12X
- Highlight: **opto edu portable colposcope, 6x portable colposcope, halogen lamp portable colposcope**



More Images



for more products please visit us on cnoec.com

Product Description

OPTO-EDU A43.1903 Portable Colposcope 6x Binocular 21V 150W Halogen Lamp

A43.1903 Colposcope

Product Description

A43.1903 Series Of Colposcope Specification						
Item	A43.1903-1	A43.1903-2	A43.1903-3			
Observation angle	Straight	45° tilted	Straight & 45° tilted (with 45° adapter)			
Magnification of binocular	6X					
Interpupillary distance	50mm-80mm					
Diopter	±5D					
Magnification changer	3-step magnification change: 0.6X,1X,1.6X					
Focal length of objective	F=200 & F=300mm (M45x0.75mm)					
Total magnification	3X, 5X, 8X, 4.7X, 7.5X, 12X					
Linear field	60.8mm, 37.9mm, 23.6mm, 40.6mm, 25.3mm, 15.8mm					
Fine focusing range	20mm					
Illumination system	Coaxial illumination with 21V150W halogen lamp cold light source fiber-optic cable, brightness adjustable, illumination>15000lx					
Filter	Green					
Fiber-optic cable	Length: 1300mm					
Floor stand	Column height adjustable fitted on the 3-jaws round base with casters					
Power supply	AC110V-AC220V					
A43.1903 Series Of Colposcope Optical Characteristics						
Focal length of objective	F=200mm			F=300mm		
Magnification on handwheel	1.6X	1X	0.6X	1.6X	1X	0.6X
Total magnification	12X	7.5X	4.7X	8X	5X	3X
Linear field (mm)	15.8	25.3	40.6	23.6	37.9	60.8
Exit pupil diameter (mm)	1.04	1.66		1.04	1.66	
Resolution (LP/mm)	67	44.5	29.7	47.2	31.5	23.6
A43.1903 Series Of Colposcope Optional Accessories						Code Number
Objective Lens	Objective F=250mm					A52.1901-25
Demonstrator	6x					A53.1901-1
Camera Adapter	C-Mount,For 1/3" Digital Camera					A55.1903-YSX
SLR Camera Adapter	For Nikon or Canon SLR camera					A55.1904-A

www.cnoec.com



A43.1903-2 Colposcope

www.cnoec.com



A43.1903-3 Colposcope



A43.1903 + Demonstrator + Beamsplitter + Camera Adapter + Digital Camera

Company Information

OPTO-EDU is a unique, long-term and excellent supplier of microscopes and educational instruments from China. We focus on the international market of microscopes and education instruments. At present our customers come from more than 80 countries, such as the United States, Britain, Russia, Canada, Germany, Denmark, Poland, Sweden, the United Arab Emirates, Qatar, Saudi Arabia, Egypt, Argentina, Korea, Thailand, Brazil and so on. OPTO-EDU as our trademark is well known in the professional field, our products ranked No.1 in Alibaba.com. We have been working tirelessly for the sale of high quality optical products and educational products made in China to the world market.

OUR WORKSHOP



OUR ADVANTAGES



As one of the most professional and most dynamic suppliers of optical and educational products in China, we have been concentrating on this field for more than 10 years. We have created the most professional product database of microscope & educational instrument, with 5000+ products and 500+ professional suppliers in system. We can provide our customers with the most comprehensive range of products in the related field, from the primary entry products to professional-level solutions, to meet the different needs of customers in the field of scientific research, teaching, industry, agriculture and various industries.



OUR SERVICE



Feedback Within 12 Hours



After Sale Service Assured




3 Years Warranty Time



Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China