



USB 3.0 CCD Microscope Hd Digital Camera 2.8MP Resolution OPTO-EDU A59.3516

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.3516
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE|Rohs
- Product Name: Digital Camera
- Resolution: 2.8MP
- Output: USB 3.0
- Sensor: CCD
- Highlight: microscope accessory, parts of the microscope, compound microscope parts



More Images



Product Description

OPTO-EDU A59.3516 USB 3.0 2.8MP CCD Microscope Digital Camera

A59.3516 Digital Camera	
Senser	SONY ICX687AQA, 2.83MP
Resolution	1940*1460@25, 1080@30fps
Pixel Size	3.69um*3.69um
Data Output	16 Bits, YUV422
Data Format	YUV422, YUYV
Partial Image	ROI
Image Processing	GAMMA, LUT, Aaturation, Acutance
Eletronic Shutter(us)	25fps:10-40000, 30fps:10-33000, step:1; manual/auto
Gain(db)	0-36, step:0.1; manual/auto
White Balance(db)	0-12, step:0.1; manual/auto
Gamma	0.5-4.00; manual/auto
Digial Interface	USB 3.0
GPIO	8pin, HIROSE HR25:PWR, 1XOPI, 1XOPO, BIOX2 ;IO support mapping
Trigger	External trigger/software trigger
Power	DC5V, <4.5W, GPIO/USB
Len Interface	C-mount
Camera Control	Drive Mode; under drive-free, low transmission
Temperature	0-50

Product Description



A59.3515 Digital Camera, USB 3.0, 16.0MP



 0086 13911110627  sale@optoedu.com  cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China