



Adjustable Educational Halogen Lamp Microscope A11.0214 With Brightness Control

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.0214
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: School
- Certification: CE|Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp|LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
binocular biological microscope,
biological compound microscope**

for more products please visit us on cnoec.com

Product Description

Educational Biological Microscope Halogen Lamp Microscopes A11.0214

Description:


This microscope is the conveniently operated optical instruments with excellent functions and new pattern, it is very suitable for usage in medical test, teaching and micro observation in hospitals and schools.


- With achromatic or semi-plan objectives
- Infinity optical system available!
- Coaxial coarse/fine focusing system with tension adjustable and up stop
- 6V 20W halogen lamp with brightness control

Specifications:

Specification	A11.0214-A	A11.0214-B
Head	Binocular head inclination 30°	
Eyepiece	WF10X(Φ18mm)	
Objective	Achromatic 4X/0.10	Semi-plan 4X/0.10
	Achromatic 10X/0.25	Semi-plan 10X/0.25
	Achromatic 40X/0.65(Spring)	Semi-plan 40X/0.65(Spring)
	Achromatic 100X/1.25(Spring, oil)	Semi-plan 100X/1.25(Spring, oil)
Focus System	Coaxial coarse/fine focus system, Minimum division of fine focusing: 2um With tensional adjustable and up stop	
Nosepiece	Quadruple Backward ball bearing inner locating	
Stage	Double layer mechanical, size 160x140mm, moving range 75x50mm	
Condenser	Abbe N.A.1.25 Rack & pinion adjustment	
Filter	Blue filter	
	Frosted filter	
Collector	For halogen lamp	
Light source	6V 20W halogen lamp with brightness control	
Optional Accessories		Item No.
Eyepiece	Wide field WF16x/11mm	A51.0203-16A
	Dividing 10x, 0.1mm/Div	A51.0205-10
Objective	Plan achromatic 4x/0.10	A52.0203-4
	Plan achromatic 5x/0.12	A52.0203-5
	Plan achromatic 10x/0.25	A52.0203-10
	Plan achromatic 20x/0.40	A52.0203-20
	Plan achromatic 40x/0.65(Spring)	A52.0203-40
	Plan achromatic 60x/0.85(Spring)	A52.0203-60
	Plan achromatic 100x/1.25(Spring, Oil)	A52.0203-100
Head	Binocular, inclination of 45°	A53.0206
	Trinocular, inclination of 45°	A53.0209
	Trinocular, inclination of 30°	A53.0207-5
	Dual viewing, inclination of 30°	A53.0203-2
Nosepiece	Quintuple (Backward ball bearing inner locating)	A54.0203-1
Filter	Green filter	A54.0209-2
	Yellow filter	A54.0209-3
Light Source	LED collector, bright and white light	A56.0204-9
Condenser	Dark field, dry	A5D.0201
	Dark field, immersion	A5D.0202
Phase Contrast Attachment	Centering telescope	A5C.0204-1
	Objective, Phase contrast plan 10x/0.25 PHP	A5C.0211-10
	Objective, Phase contrast plan 20x/0.40 PHP(Spring)	A5C.0211-20
	Objective, Phase contrast plan 40x/0.65 PHP(Spring)	A5C.0211-40
	Objective, Phase contrast plan 100x/1.25 PHP (Spring,oil)	A5C.0211-100
	PH-I Turnplate phase contrast condenser	A5C.0202-1

	PH-II Turnplate phase contrast condenser	A5C.0202-4
	Flashboard phase contrast set	A5C.0202-2
	Pull board phase contrast set	A5C.0202-3
Infinity Optical System	Binocular head inclination of 30°	A53.0205-2
	Trinocular head inclination of 30°	A53.0208-2
	Centering telescope	A5C.0204-2
	Objective, infinity plan achromatic 4x/0.10	A52.0206-4
	Objective, infinity plan achromatic 10x/0.25	A52.0206-10
	Objective, infinity plan achromatic 20x/0.25	A52.0206-20
	Objective, infinity plan achromatic 40x/0.65(Spring)	A52.0206-40
	Objective, infinity plan achromatic 100x/1.25(Spring, oil)	A52.0206-100
	Objective, infinity phase contrast plan 10x/0.25 PHP	5C.0216-10
	Objective, infinity phase contrast plan 20x/0.40 PHP (Spring)	5C.0216-20
	Objective, infinity phase contrast plan 40x/0.65 PHP (Spring)	5C.0216-40
	Objective, infinity phase contrast plan 100x/1.25 PHP (Spring, oil)	5C.0216-100
CCD adapter	0.4x	A55.0202-1
	0.5x	A55.0202-4
	1x	A55.0202-2
	0.5x with dividing 0.1mm/Div	A55.0202-3
Digital Camera Adapter	Canon digital camera adapter (A610, A620, A630, A640)	A55.0204-01
Photo Adapter	2.5x/4x change over photograph attachment with 10x viewing eyepiece	A55.0201-1
	4x focusing photograph attachment	A55.0201-2
	MD adapter	A55.0201-3
	PK adapter	A55.0201-4
Adapter	Input: 230V<50/60HZ Output: 9V 500mA	A56.0205-1
	Input: 110V<50/60HZ Output: 9V 500mA	A56.0205-2
Polarizing Attachment	Polarizer 360° rotatable	A5P.0201
	Analyzer	





Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China