



A11.1009 1000X LED Light Source Biological Microscope For Student / School

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1009
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student|school
- Certification: CE|Rohs
- Drawtube: Binocular|Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
binocular biological microscope,
biological compound microscope**



for more products please visit us on cnoec.com

Product Description

A11.1009 1000X LED light source Biological Microscope

High Quality Compact Design Student Biological Microscope;
WF10x/18mm Wide Field Eyepiece;
Monocular, Binocular Head 45° Inclined, 360° Rotatable;
Double Layer Mechanical Stage 110x120mm Moving 70x20mm;
Coaxial Coarse & Fine Focusing Konb, Light Source Brightness Adjustable.



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085


A11.1009

Student Biological Microscope





A11.1009-D Monocular Head



A11.1009-E Binocular Head




Achromatic Objectives



Double Layer Mechanical Stage



Brightness Adjustable



Abbe Condenser N.A.1.2, With Iris Diaphragm



A11.1009-D



A11.1009-E

| Items | A11.1009 Student Biological Microscope | M | B |
|----------------------|---|---|---|
| Head | Monocular Head, 45° Inclined, 360° Rotatable | ● | |
| | Sliding Binocular Head, 45° Inclined, 360° Rotatable | | ● |
| Eyepiece | Wide Field WF10x/18mm | ● | ● |
| | WF5x | ○ | ○ |
| | WF16x | ○ | ○ |
| | WF10x Eyepiece With Pointer | ○ | ○ |
| Nosepiece | Triple Nosepiece | ● | |
| | Quadruple Nosepiece | | ● |
| Achromatic Objective | 4x/0.10, 10x/0.25, 40x/0.65(S) | ● | |
| | 4x/0.10, 10x/0.25, 40x/0.65(S), 100x/1.25(S,Oil) | | ● |
| | 20x 60x | ○ | ○ |
| Magnification | 40x~400x | ● | |
| | 40x~1000x | | ● |
| Stage | Double Layer Mechanical Stage, 110*120mm, Moving Range 70*20mm | ● | ● |
| Focusing | Coaxial Coarse & Fine Focusing | ● | ● |
| Condenser | Abbe Condenser N.A.1.2, With Iris Diaphragm, With Blue Filter, Rack & Pinion Up/Down Adjustment | ● | ● |
| Light Source | Incandescent Lamp Illumination 15W, Brightness Adjustable | ● | |
| | LED Lamp Illumination 1W, Brightness Adjustable | | ● |
| | Mirror | ○ | ○ |
| Spare Parts | Dust Cover, Spare Bulb & Fuse | ● | ● |

Note: "●" In Table Is Standard Outfits, "○" Is Optional Accessories.



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China