



## 1000X Binocular Handheld Digital Microscope OPTO-EDU A11.1133 CE / Rohs Certification

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1133
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Student|School|Hobby|Home
- Certification: CE|Rohs
- Drawtube: Monocular/Binocular
- Illumination: LED
- Max Magnification: 1000X
- Theory: Biological Microscope
- Highlight: **high powered microscope ,  
binocular biological microscope,  
biological compound microscope**



### More Images



## Product Description

### Portable Monocular/Binocular Biological Microscope With Coaxial Coarse And Fine Focusing A11.1133

#### Specifications:

A11.1133 Biological Microscope Specification		A	B	C	D
Head	Monocular head, 30° inclined, 360° rotatable	●	●		
	Sliding binocular head, 45° inclined, 360° rotatable			●	
	Seidentopf binocular head, 30° inclined, 360° rotatable				●
Eyepiece	WF10X with pointer	●	●	●	●
Objective	Achromatic 4X	●	●	●	●
	Achromatic 10X	●	●	●	●
	Achromatic 40X, spring	●	●	●	●
	Achromatic 100X, spring/oil	●	●	●	●
Nosepiece	Triple revolving nosepiece	●			
	Quadruple revolving nosepiece		●	●	●
Stage	Plain stage with clips, size 100x90mm	●			
	Attachable mechanical stage, size 100x90mm		●	●	
	Double layer Mechanical stage, size 100x90mm				●
Condenser	Single lens condenser with Disc diaphragm	●			
	N.A. 1.25 condenser with Iris diaphragm		●	●	●
Focusing	Coaxial coarse and fine focusing knobs	●	●	●	●
Light Source	LED lampe with rechargeable battery	●	●	●	●
	Power supply adapter	●	●	●	●





**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China