



40x - 400x Biological Monocular Optical Microscope A11.1506-A1 Simple Structure

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1506
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Total Magnification: 40x-400x
- Head: Monocular Vertical Head Tube
- Eyepiece: WF10x
- Objectives: 4x, 10x, 40x
- Focusing: Coarse Adjustment Range: 20mm
- Stage: Plain Stage With Clips, 100mmx90mm
- Illuminator: Concave-Plano Mirror Dia. 40mm
- Highlight: **high powered microscope, binocular biological microscope, biological compound microscope**



More Images



Product Description

OPTO-EDU A11.1506-A1 Monocular Cheapest professional educational student microscope

Total Magnification: 40X~400X
 Monocular Vertical Tube, Coarse Focus
 Simple Structure And Nice Appearance
 Featured By Clear Image, Easy Operation


The Biological Microscope Model A11.1506, Can Be Used In Teaching Experiment For Primary And Secondary Schools. It Is Featured By Clear Image, Easy Operation, Simple But Reliable Structure And Nice Appearance. Stration And Other Occasions.

A11.1506 Specification		A1	A2	A3
Total Magnification	40x-400x	●	●	●
Head	Monocular Vertical Head Tube	●	●	-
	Monocular 45°Inclined,360°Rotatable	-	-	●
Eyepiece	WF10x	●	●	●
Objectives	4x	●	●	●
	10x	●	●	●
	40xs	●	●	●
Focusing	Coarse Adjustment Range: 20mm	●	-	-
	Seperate Coarse And Fine Adjustment,Coarse Adjustment Range: 20mm, Fine Adjustment Range : 1.4~1.6mm	-	●	●
Stage	Plain Stage With Clips,100mmx90mm	●	●	●
Illuminator	Concave-Plano Mirror Dia. 40mm	●	●	●
Optional Accessories				
Illuminator	Electrical Illuminator (110V/220V 20W)			

A11.1506 BIOLOGICAL MICROSCOPE



A11.1506 Biological Microscope



PACKING AND SHIPPING

Packing : Cardbox And Polyfoam + Strong Master Carton 15*21*34cm, G.W.: 2.0Kgs
Shipping Time : 3-10 Days From Stock



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China