



## WF10x 16x Up Bottom LED Light Biological Microscope A11.1545 For Student

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1545
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Head: Monocular Head, 45° Inclined
- Eyepiece: WF10x,16x
- Objectives: 4x,10x,40xs
- Focusing: Coarse Focusing Knob, Focusing Range: 20mm
- Stage: 90mm\*90mm, Clip Stage
- Condenser: Single Lens Condenser, 6 Holes Disc Diaphragm
- Light Source: Up & Bottom LED Light, Both Brightness Adjustable
- Highlight: **high powered microscope, binocular biological microscope, biological compound microscope**



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

WF10x,16x OPTO-EDU Up & Bottom LED Light biological Monocular Student Microscope A11.1545

A11.1545 Biological Microscope	
Specification	A11.1545 A11.1545-E
Digital eyepiece	0.3M
Head	Monocular Head, 45° Inclined
Eyepiece	WF10x,16x
Objectives	4x,10x,40xs
Focusing	Coarse Focusing Knob, Focusing Range: 20mm
Stage	90mm*90mm, Clip Stage
Condenser	Single Lens Condenser, 6 Holes Disc Diaphragm
Light Source	Up & Bottom LED Light, Both Brightness Adjustable
Power	Cordless Batteries (3xAAA Batteries NOT Included) Adapter Included (110V/220V Optional)
Standard Accessories	
Packing Box	Packed In Suitcase Packing Box
Experimentation Appliance	Plastic Microtom 1 pc
	Hatching Box 1 pc
	Tweezer 1pc
	Sea Salt 1 Bottle
	Yeast 1 Bottle
	Worm Eggs 1 Bottle
	Glue 1 Bottle
Blank Slides 5 pcs, Prepared Slides 5 pcs, In 1 Box	

A11.1545 Biological Microscope	
Specification	A11.1545 A11.1545-E
Digital eyepiece	0.3M
Head	Monocular Head, 45° Inclined
Eyepiece	WF10x,16x
Objectives	4x,10x,40xs
Focusing	Coarse Focusing Knob, Focusing Range: 20mm
Stage	90mm*90mm, Clip Stage
Condenser	Single Lens Condenser, 6 Holes Disc Diaphragm
Light Source	Up & Bottom LED Light, Both Brightness Adjustable
Power	Cordless Batteries (3xAAA Batteries NOT Included) Adapter Included (110V/220V Optional)
Standard Accessories	
Packing Box	Packed In Suitcase Packing Box
Experimentation Appliance	Plastic Microtom 1 pc
	Hatching Box 1 pc
	Tweezer 1pc
	Sea Salt 1 Bottle
	Yeast 1 Bottle
	Worm Eggs 1 Bottle
	Glue 1 Bottle
Blank Slides 5 pcs, Prepared Slides 5 pcs, In 1 Box	



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China