



## Binocular 45° Inclined Head Portable Digital Microscope 3.36x - 216x Total Magnification

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1826
- Minimum Order Quantity: 1 pc
- Price: Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Total Magnification: 3.36x-216x
- Head: Binocular 45° Inclined, 48-75mm
- Eyepiece: WF10X/22mm, WF20X/12mm
- Zoom Lens: 0.7-4.5x
- Bridge Magnification: 1.2x
- Auxiliary Lens: 0.4x, 2.0x
- Highlight: **bullet comparison microscope, binocular optical microscope**



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### OPTO-EDU A18.1826 3.36~216x Binocular 45° Inclined Forensic bullet Comparison Microscope

A18.1826 Comparison Microscope	
T o t a l M a g n i f i c a t i o n	3.36x-216x
H e a d	Binocular 45° Inclined, 48-75mm
E y e p i e c e	WF10X/22mm, WF20X/12mm
Z o o m L e n s	0.7-4.5x
B r i d g e M a g n i f i c a t i o n	1.2x
A u x i l i a r y L e n s	0.4x, 2.0x
V i e w F i e l d	Full Left, Full Right, Split Image, Overlapping Image, 2mm-60mm
W o r k i n g	Moving Range: X-54mm, Y-54mm, Z-54mm, 0°-45° Inclination Adjustable
	Joint Left/Right Moving Range 62mm, Joint Up/Down Range 62mm


Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

Stage	Round Stage, 360° Rotating
	Bullet Holder 2 Sets
Light Source	5W LED+12V/50W Halogen
	Coaxial Light
	Polarizing Light
	UV Filter
	Filter Red/Blue/Yellow/Green
Adapter	0.5x C-Mount
	Photo Adapter For Cannon SLR
Function Panel	
	Touch Type, User-Friendly Design, Easy To Operate And Accurate
Power	Power Supply 80-240V

**OPTO-EDU**

OPTO-EDU (BEIJING) CO.,LTD.



## A18.1826 Comparison Microscope

### HD Image

Visual HD images, static data acquisition via camera installation

### Application

It is the main scientific and technical means of national criminal investigation technology, procuratorial technology and forensic identification, and has been widely used in the field of judicial identification technology.



0.4x Auxiliary Lens



Bullet Holder



Filter Red/Blue/Yellow/Green/UV

## Product Real Shot

### ■ Coin Test Switching Field Of View



Split Image



Overlapping Image



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China