



## 0.7 - 4.5x Objective Lens Compound Optical Microscope Cylindrical Trace Extending Camera A18.1837

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1837
- Minimum Order Quantity: 1 pc
- Price: Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Total Magnification: 0.14~50x
- Objective Lens: 0.7~4.5x
- Focusing Range: 300mm
- Detection Range: Dia. 3~120 Mm
- Height: 2~300 Mm
- Unfolding Rngle: 0~360 Degree Customizable
- Light Source: 5W LED
- Camera: Internal & External Dual Cameras
- Software: Including Software+Laptop
- Highlight: **bullet comparison microscope, binocular optical microscope**



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

OPTO-EDU A18.1837 0.14~50x Objective lens 0.7~4.5x Cylindrical Trace Extending Camera

A18.1837 Cylindrical Trace Extending Camera	
Total Magnification	0.14~50x
Objective lens	0.7~4.5x
Focusing Range	300mm
Detection Range	Dia. 3~120 mm
Height	2~300 mm
Unfolding Rngle	0~360 Degree Customizable
Light Source	5W LED
Camera	Internal & External Dual Cameras
Software	Including Software+Laptop



Our Product Introduction


for more products please visit us on [cnoec.com](http://cnoec.com)



**Opto-Edu (Beijing) Co., Ltd.**

 0086 13911110627

 [sale@optoedu.com](mailto:sale@optoedu.com)

 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China