



A18.1838 0.14 - 50x Forensic Comparison Microscope Cylindrical Trace Extending Camera

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1838
- Minimum Order Quantity: 1 pc
- Price: Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Total Magnification: 0.14~50x
- Objective Lens: 0.7~4.5x
- Focusing Range: 300mm
- Detection Range: Dia. 3~120 Mm
- Height: 2~300 Mm
- Unfolding Rngle: 0~360 Degree Customizable
- Light Source: Multiband Diffuse Light Source
- Camera: Internal & External Dual Cameras
- Software: Including Software+Laptop
- Highlight: **bullet comparison microscope, binocular optical microscope**



Product Description

OPTO-EDU A18.1838 0.14~50x Detection Range Dia. 3~120 mm Cylindrical Trace Extending Camera

A18.1838 Cylindrical Trace Extending Camera	
Total Magnification	0.14~50x
Objective lens	0.7~4.5x
Focusing Range	300mm
Detection Range	Dia. 3~120 mm
Height	2~300 mm
Unfolding Rngle	0~360 Degree Customizable
Light Source	Multiband Diffuse Light Source
Camera	Internal & External Dual Cameras
Software	Including Software+Laptop

--	--	--	--	--	--	--

Our Product Introduction

for more products please visit us on cnoec.com



Opto-Edu (Beijing) Co., Ltd.

0086 13911110627

sale@optoedu.com

cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China