



Two Camera Rings Multi Spectral Observation Camera OPTO-EDU A18.1843 High Precision

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1843
- Minimum Order Quantity: 1 pc
- Price: Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Spectrum1: White Light
- Spectrum2: Purple Light 415nm
- Spectrum3: Blue Light 450nm
- Spectrum4: Green Light 530nm
- Spectrum5: Yellow Light 580nm
- Spectrum6: Red Light 615nm
- Camera: Two Camera Rings
- Highlight: **bullet comparison microscope, binocular optical microscope**



for more products please visit us on cnoec.com

Product Description

OPTO-EDU A18.1843 Two Camera Rings Multispectral Observation Camera

A18.1843 Multispectral Observation Camera	
S p e c i f i c a t i o n s	White Light
	Purple Light 415nm
	Blue Light 450nm
	Green Light 530nm
	Yellow Light 580nm
	Red Light 615nm
C a m e r a	Two Camera Rings



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China

Our Product Introduction

for more products please visit us on cnoec.com