



LCD Standalone Inspection Digital Microscope 600x Brightness Adjustable A33.5006

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A33.5006
- Minimum Order Quantity: 1 pc
- Price: Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month




Product Specification

- Magnification Ratio: 600x
- Still Capture Resolution: 8M, 5M, 3M, 1.3M, VGA
- Focus Range: Manual Focus From 5mm To 80mm
- Frame Rate: Max 30f/s Under 600 Lus Brightness
- Video Format: AVI
- Photo Format: JPEG
- Light Source: LEDs (Brightness Adjustable)
- Highlight: **digital trinocular microscope,
usb digital microscope,
infinity plan microscope**



More Images




Product Description

OPTO-EDU Digital Microscope 4.3-inch LCD 600x LED Brightness Adjustable A33.5006

The 4.3-inch LCD Standalone Inspection Digital Microscope is a fantastic microscope that provides a good way for you to share high definition microscopic observation with others. The 4.3" quality LCD panel brings vivid crystal image. Equipped with professional microscopic lens, it helps you to see tiny objects easily. The magnification can be up to 600x. Built-in rechargeable Li-ion battery provides a 6-7rs standalone working time. The 8 bright LEDs, with brightness adjustable, are bright enough for perfect observation.

A33.5006 4.3-inch LCD Inspection Microscope Specification	
Magnification Ratio	600x
Still Capture Resolution	8M, 5M, 3M, 1.3M, VGA
Focus Range	Manual Focus From 5mm To 80mm
Frame Rate	Max 30f/s Under 600 Lus Brightness
Video Format	AVI
Photo Format	JPEG
Light Source	LEDs (Brightness Adjustable)
Battery	Built-in Li-ion Battery
DC Input	5V/1A
TV-out	To Any Monitor With TV-in
Storage	Micro-SD Card
Bundle Software	Portable Capture Plus With Measurement
Operation System For Bundle Software	Windows 7/8/10 Mac OS X 10.12 Or Above Android Smartphone/Tablet
Product Size	116mm (L) x 95mm (W)
Weight	Microscope-220g; Stand-220g
Package Contents	Microscope Observation Stand Software CD USB Cable TV Cable Micro-SD Card (1GB) User's Manual

PART NAME AND FUNCTION



1	Focus Knob	9	Micro-SD Card Slot
2	Menu	10	Digital Transform Multiples
3	Video, Photo Capture, Playback Switch	11	Battery Charging Indicator Light (Red For Charging, Full Of Rear Light Out)
4	Up	12	Television Outlet
5	Down	13	Charging Port: Charge The Internal Battery USB Connected Computer: Read The Card, Or Enable Live Video On Your PC/Android Device
6	Video & Photo Capture		
7	Boot to Turn It Off		
8	LED Source Adjust		



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China