



## Eyepiece WF10x Stereo Optical Microscope A23.1122 LED Illumination Zoom Lens 1 - 4.5x

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.1122
- Minimum Order Quantity: 1 pc
- Price: Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Eyepiece: WF10x
- Objective: Zoom Lens 1~4.5x
- Stand: With Upper/Bottom Light Source
- Incident Illumination: 3W LED Brightness Adjustable
- Transmitted Illumination: 1W LED Brightness Adjustable
- Head: Binocular & Trinocular
- Highlight: **stereo zoom microscopes,  
stereo binocular microscopes,  
stereo dissecting microscope**



### More Images



## Product Description

### Eyepiece WF10x Stereo Optical Microscope A23.1122 LED Illumination Zoom Lens 1 - 4.5x

A23.1122 Stereo Microscope, 1.0~4.5x, 1:4.5		B	T	AB	AT
Head	Binocular, Interpupillary Distance 55~75mm	•		•	
	Trinocular, Interpupillary Distance 55~75mm		•		•
Eyepiece	WF10x/22mm	•	•	•	•
Objective	Zoom Lens 1~4.5x, Zoom Ratio:1:4.5, Working Distance120mm	•	•	•	•
Stand	Track Stand With Handle, With Upper/Bottom Light Source	•	•		
	Pole Stand With Handle, With Upper/Bottom Light Source			•	•
Illumination	Incident Illumination 3W LED Brightness Adjustable	•	•	•	•
	Transmitted Illumination 1W LED Brightness Adjustable	•	•	•	•
Optional Accessory	Eyepiece: WF10/23mm,WF15/17mm,WF20/12.5mm				
	Right Light Source				
	Auxiliary Objective 0.5x, 2x				
	C-mount				


THE KEY TO MICRO WORLD

OPTO-EDU



# A23.1122

Stereo Microscope, 1.0~4.5x, 1:4.5



Binocular, Trinocular  
Head With  
WF10x/22mm

Up LED 1W



BottomLED 3W



## Product Specification



A23.1122-B



A23.1122-AB



A23.1122-T



A23.1122-AT



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China