



LED Digital Portable Metallurgical Optical Microscope A31.2501 3.1M 100 - 400x

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A31.2501
- Minimum Order Quantity: 1 pc
- Price: Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Total Magnification: 100-400x
- Mechanical Tube Length: 160mm
- Head: Monocular Vertical Head
- Eyepiece: WF10X
- Focusing: Coaxial Coarse & Fine Focus Knob, Focusing Range 36mm
- Light Source: LED Coaxial Illumination, Brightness Adjustable
- Highlight: **reflected light microscope, halogen lamp microscope**


A31.2501 Digital Portable Metallurgical Microscope



A50 2207-3.1M
3.1M Digital Camera



Compact All-in-One
Portable Design For
Checking Metal
Surface In Anywhere



Supply In Strong
Aluminum Case, Magnet
Stand Optional



More Images



Product Description

100 - 500x Digital Portable Metallurgical Optical Microscope A31.2501

Small & Simple, with bright white LED light built-in, and high quality image!

| | | |
|---|--|--------------|
| A31.2501-B Digital Portable Metallurgical Microscope | | |
| Total Magnification | 100-400x | |
| Mechanical Tube Length | 160mm | |
| Head | Monocular Vertical Head | |
| Eyepiece | WF10X | |
| LWD Plan Objective (No Cover Glass) | PLL 10X/0.25 W.D. 7.3mm | |
| | PLL 40X/0.65 W.D. 0.5mm | |
| Focusing | Coaxial Coarse & Fine Focus Knob, Focusing Range 36mm | |
| Light Source | LED Coaxial Illumination, Brightness Adjustable, Powered By Rechargeable Batteries (Battery NOT Included For Air Shipment) | |
| Packing | Portable Aluminium Case | |
| Digital Camera | 3.1M/MT9T001(C) 1/2" (6.55x4.92) Max.Resolution 2048/*1536 Pixel Size 3.2 x 3.2µm Responsivity 1.0 V/lux-sec 2048*1536/10fps 1024*768/22fps 680*510/43fps Dynamic Range 61dB SN Ration 43dB Binning 1x1,2x2,3x3 Exposure 0.244~2000ms | |
| A31.2501-B Digital Portable Metallurgical Microscope Optional Accessories | | |
| Adapter | Photo Adapter For SLR Digital Camera Canon Series | A55.2622-C |
| | Photo Adapter For SLR Digital Camera Nikon Series | A55.2622-N |
| Polarizing Set | Polarizer & Analysiser Slides | A5P.2601 |
| LWD Plan Objective (No Cover Glass) | 20x/0.40mm, WD8.6 | A5M.0234-20 |
| | 50x/0.70mm, WD2.20, Spring | A5M.0234-50A |
| | 60x/0.75mm, WD1.34, Spring | A5M.0234-60 |
| | 80x/0.80mm, WD0.96, Spring | A5M.0234-80A |
| | 100x/1.25mm, WD0.33, Spring, Oil | A5M.0234-100 |
| Magnet Stand | Magnet Stand, X/Y Moving Range 10x10mm | A54.2501 |

THE KEY TO MICRO WORLD

OPTO-EDU





A59.2207-3.1M CCD Camera



Portable Aluminium Case



Total Magnification 100-400x




No Cover Glass LWD Plan Objective (Optional)



LED Coaxial Illumination, Brightness Adjustable



Coaxial Coarse & Fine Focus



Polarizer & Analyser Slides (Optional)




Rechargeable Power




Photo Adapter For SLR Digital Camera (Optional)





Magnet Stand (Optional)



www.optoedu.com



4




sale@cnoec.com



Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China