



## LED OPTO-EDU A11.2601-BP Laboratory Microscope Binocular

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.2601-BP
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Head: Seidentopf Binocular Head Inclined 30°
- Applications: Student|school
- Certification: CE|Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp|LED
- Max Magnification: 1000X
- Product Name: Microscope
- Highlight: **OPTO-EDU laboratory microscope, A11.2601-BP laboratory microscope, OPTO-EDU microscope lab**



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

A11.2601 Biological Microscope is designed for teaching clinical verify. This series microscope with original style, steady structure, convenient operation and clear image, is suitable for observing various biological specimens. And mostly applied in school, college, hospital and family.



### Biological Microscope CNOPTEC Series

#### A11.2601

- ◆ Wide Field Eyepiece WF10x/18mm
- ◆ High Contrast Achromatic Objective 4x, 10x, 40x, 100x
- ◆ Mechanical Stage 115x125mm, Moving Range 70x35mm
- ◆ Full Metal Body High Rigid Construction
- ◆ 3W LED, Cordless Rechargeable Function Optional



Pretty

Durable

Good Effect

A11.2601 Biological Microscope Is Designed For Teaching, Clinical Verify. The Series Microscope With Original Style, Steady Structure, Convenient Operation And Clear Image Is Suitable For Observing Various Biological Specimens, And Mostly Applied In School, College, Hospital And Family.

A11.2601 Biological Microscope is designed for teaching clinical verify. This series microscope with original style, steady structure, convenient operation and clear image, is suitable for observing various biological specimens. And mostly applied in school, college, hospital and family.



## **01 Head**

WF10x-18mm  
Eyepiece Is Anti-fungal, Anti-fungal System  
Built In Head

## **02 Objective**

Achromatic Objective: 4x, 10x, 40x(s),  
100x/1.25(oil)(s)  
Objective Is Anti-fungal



## **03 Stage**

Mechanical Stage  
Stage Size:125mm×115mm,  
Travel: 35mm×75mm

## **04 Focus**

Coaxial Coarse & Fine Focus  
Stage Moving Knob Can Move Stage  
Turn Left Or Right



## **05 Condenser**

Abbe Condenser  
N.A. 1.25, Adjustable Aperture, Aperture  
Center Is Adjustable  
Built In Blue Filter,  $\phi$ 38.5mm

## **06 Light Source**

LED Illumination Systems (Non-recharge-  
able)  
Rotation Potentiometer With Power Switch





Items	A11.2601 Student Biological Microscope	-M	-B	-T	-BP	-TP	Cata. No.
Head	Monocular Head Inclined 30°, Rotatable 360°	●					A53.2601-M
	Seidentopf Binocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm		●		●		A53.2601-B
	Seidentopf Trinocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm, Light Split Fixed P20/E80			●		●	A53.2601-T
Eyepiece	WF10x/18mm	●	●●	●●	●●	●●	A51.2602-1018
	WF16x/13mm	○	○	○	○	○	A51.2602-1613
	WF10x/18mm With Pionter	○	○	○	○	○	A51.2604-1018
	WF10x/18mm, With Micrometer X-0.1mm*180 Div / Y-Single Line	○	○	○	○	○	A51.2605-1018
	WF10x/18mm With Micrometer X-0.1mm*180 Div / Y-Single Line, Diopter Adjustable	○	○	○	○	○	A51.2605-1018D
Nosepiece	Quadruple, Backward	●	●	●	●	●	A54.2610-N03
Achromatic Objective	4x/0.1, W.D.=18.6mm	●	●	●			A52.2601-4
	10x/0.25, W.D.=6.5mm	●	●	●			A52.2601-10
	20x/0.40, W.D.=1.75mm	○	○	○			A52.2601-20
	40x/0.6(S), W.D.=0.47mm	●	●	●			A52.2601-40
	60x/0.80(S), W.D.=0.10mm	○	○	○			A52.2601-60
	100x/1.25(S,O), W.D.=0.07mm	●	●	●			A52.2601-100
Plan Objective	4x/0.1, W.D.=14.5mm				●	●	A52.2603-4
	10x/0.25, W.D.=5.65mm				●	●	A52.2603-10
	20x/0.40, W.D.=5.65mm				○	○	A52.2603-20
	40x/0.6(S), W.D.=0.85mm				●	●	A52.2603-40
	100x/1.25(S,O), W.D.=0.07mm				●	●	A52.2603-100
Stage	Double Layer Mechanical Size 125×115mm, Travel Range 35×75mm	●	●	●	●	●	A54.2601-S03
Condenser	Abbe N.A. 1.25(Iris Diaphragm)	●	●	●	●	●	A56.2614-09
Focusing	Coaxial Coarse And Fine Focusing Knobs, Travel Range: 22mm, Scale: 2μm	●	●	●	●	●	
C-mount	1.0x (Focus Adjustable)			○		○	A55.2601-B10
	0.5x (Focus Adjustable)			○		○	A55.2601-B05
Light Source	3W LED Illumination Systems	●	●	●	●	●	A56.2613-3W
	3W LED Illumination Systems (Rechargeable)	○	○	○	○	○	A56.2651-3WR

	12V 20W Halogen Lamp, Wide Voltage: 100V-240V	○	○	○	○	○	A56.2650-20W
	Plano-Concave Mirror	○	○	○	○	○	A56.2601-B
Dark Field	Dark Field Ring, For 4x~40x Objective	○	○	○	○	○	A5D.2611-B
Polarizing	Analyzer + Polarizer Set	○	○	○	○	○	A5P.2601-B
Package	Carton Size 337*242*417mm, 1 Unit/Ctn						
Weight	Gross Weight 7.0kgs, Net Weight 5.5kgs						
Note: "●" In Table Is Standard outfits, "○" Is Optional Accessories.							



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China