



1000X OPTO-EDU A11.2601-M Educational Microscope

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.2601-M
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Head: Monocular Head Inclined 30°, Rotatable 360°
- Applications: Student|school
- Certification: CE|Rohs
- Drawtube: Monocular
- Illumination: Halogen Lamp|LED
- Max Magnification: 1000X
- Product Name: Microscope
- Highlight: **OPTO-EDU Educational Microscope, 1000X Educational Microscope, 1000X Educational Biological Microscope**



More Images



for more products please visit us on cnoec.com

Product Description

1000X Elementary Biological Microscope Research , Rohs CE A11.2601

Description:

A11.2601 Student Microscope, modern design with high quality optical system, supplying clear image, is good to be used in high school, university for student.

Wide Field Eyepiece WF10x/18mm



Biological Microscope CNOPTEC Series

A11.2601

- ◆ Wide Field Eyepiece WF10x/18mm
- ◆ High Contrast Achromatic Objective 4x,10x,40x,100x
- ◆ Mechanical Stage 115x125mm, Moving Range 70x35mm
- ◆ Full Metal Body High Rigid Construction
- ◆ 3W LED, Cordless Rechargeable Function Optional



Pretty

Durable

Good Effect

A11.2601 Biological Microscope Is Designed For Teaching, Clinical Verify. The Series Microscope With Original Style, Steady Structure, Convenient Operation And Clear Image Is Suitable For Observing Various Biological Specimens, And Mostly Applied In School, College, Hospital And Family.



01 Head

WF10x-18mm
Eyepiece Is Anti-fungal, Anti-fungal System
Built In Head

02 Objective

Achromatic Objective: 4x, 10x, 40x(s),
100x/1.25(oil)(s)
Objective Is Anti-fungal



03 Stage

Mechanical Stage
Stage Size:125mm×115mm,
Travel: 35mm×75mm

04 Focus

Coaxial Coarse & Fine Focus
Stage Moving Knob Can Move Stage
Turn Left Or Right



05 Condenser

Abbe Condenser
N.A. 1.25, Adjustable Aperture, Aperture
Center Is Adjustable
Built In Blue Filter, ϕ 38.5mm

06 Light Source

LED Illumination Systems (Non-recharge-
able)
Rotation Potentiometer With Power Switch





Items	A11.2601 Student Biological Microscope	-M	-B	-T	-BP	-TP	Cata. No.
Head	Monocular Head Inclined 30°, Rotatable 360°	●					A53.2601-M
	Seidentopf Binocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm		●		●		A53.2601-B
	Seidentopf Trinocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm, Light Split Fixed P20/E80			●		●	A53.2601-T
Eyepiece	WF10x/18mm	●	●●	●●	●●	●●	A51.2602-1018
	WF16x/13mm	○	○	○	○	○	A51.2602-1613
	WF10x/18mm With Pionter	○	○	○	○	○	A51.2604-1018
	WF10x/18mm, With Micrometer X-0.1mm*180 Div / Y-Single Line	○	○	○	○	○	A51.2605-1018
	WF10x/18mm With Micrometer X-0.1mm*180 Div / Y-Single Line, Diopter Adjustable	○	○	○	○	○	A51.2605-1018D
Nosepiece	Quadruple, Backward	●	●	●	●	●	A54.2610-N03
Achromatic Objective	4x/0.1, W.D.=18.6mm	●	●	●			A52.2601-4
	10x/0.25, W.D.=6.5mm	●	●	●			A52.2601-10
	20x/0.40, W.D.=1.75mm	○	○	○			A52.2601-20
	40x/0.6(S), W.D.=0.47mm	●	●	●			A52.2601-40
	60x/0.80(S), W.D.=0.10mm	○	○	○			A52.2601-60
	100x/1.25(S,O), W.D.=0.07mm	●	●	●			A52.2601-100
Plan Objective	4x/0.1, W.D.=14.5mm				●	●	A52.2603-4
	10x/0.25, W.D.=5.65mm				●	●	A52.2603-10
	20x/0.40, W.D.=5.65mm				○	○	A52.2603-20
	40x/0.6(S), W.D.=0.85mm				●	●	A52.2603-40
	100x/1.25(S,O), W.D.=0.07mm				●	●	A52.2603-100
Stage	Double Layer Mechanical Size 125×115mm, Travel Range 35×75mm	●	●	●	●	●	A54.2601-S03
Condenser	Abbe N.A. 1.25(Iris Diaphragm)	●	●	●	●	●	A56.2614-09
Focusing	Coaxial Coarse And Fine Focusing Knobs, Travel Range: 22mm, Scale: 2μm	●	●	●	●	●	
C-mount	1.0x (Focus Adjustable)			○		○	A55.2601-B10
	0.5x (Focus Adjustable)			○		○	A55.2601-B05
Light Source	3W LED Illumination Systems	●	●	●	●	●	A56.2613-3W
	3W LED Illumination Systems (Rechargeable)	○	○	○	○	○	A56.2651-3WR

	12V 20W Halogen Lamp, Wide Voltage: 100V-240V	○	○	○	○	○	A56.2650-20W
	Plano-Concave Mirror	○	○	○	○	○	A56.2601-B
Dark Field	Dark Field Ring, For 4x~40x Objective	○	○	○	○	○	A5D.2611-B
Polarizing	Analyzer + Polarizer Set	○	○	○	○	○	A5P.2601-B
Package	Carton Size 337*242*417mm, 1 Unit/Ctn						
Weight	Gross Weight 7.0kgs, Net Weight 5.5kgs						
Note: "●" In Table Is Standard outfits, "○" Is Optional Accessories.							



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China