



1000X 5.0M Resolution 9.7" Portable Lcd Digital Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A33.1502
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Authenticate: CE
- Export Area: Europe & North America
- Monthly Sales: 5000 Units
- Application Field: Efficient And Laboratory
- Magnification: 1000x
- Resolution: 5.0M
- Highlight: **9.7" Portable Lcd Digital Microscope, 1000X Portable Lcd Digital Microscope, 5.0M Portable Lcd Digital Microscope**



More Images



Product Description

1. 9.7" LCD, 5.0M, Trinocular, WF10x/20, Quaduple, Achromatic HC 4x,10x,40x100x, Double Layer Mechanical Graphite 143*132mm, Moving 75x50mm, Abbe N.A.1.25, Coaxial, 3W LED, 1.0x C Mount
2. 1.9.7" LCD Touch Screen Pad, With Android System + Measuring Software
3. High Eyepoint Wide Field Eyepiece WF10x/20mm Dia.30mm
4. Achromatic HC 4x,10x,40x100x Objective, Upgradeable To Infinity
5. Double Layer Mechanical Stage 143*132mm, Upgradeable To 216x150mm
6. 3W LED Kohler Illumination With Removable Light House For Easy Maintenance

OPTO-EDU

A33.1502 Research Level Digital Microscope

9.7" LCD | Android Pad | 5.0M



Our Product Introduction

for more products please visit us on cnoec.com

9.7" Android Pad + 5.0M Camera

Sharp Image | All-in-One | Photo & Video



Explore Micro World Together

Multipul Usage | Easy To Share | Nice For Teaching



Laboratory



School




Hospital





| | |
|---|---|
| A33.1502 9.7" LCD Digital Biological Microscope | |
| Head | Trinocular Head, 30° Inclined, Interpupillary Distance 48-75mm, Both Eyepiece Tube Diopter Adjustable |
| Eyepiece | WF10x/20mm, Diopter Adjustable, Dia. 30mm |
| Objective | High Contrast Chromatic Free 4x, 10x, 40x(S), 100x(Oil, S) |
| Focus Knob | Coaxial Coarse & Fine Focusing, Fine Division 0.002mm, Coarse Stroke 37.7mm/Rotation, Focusing Range 24mm |
| Nosepiece | Quadruple, Backward With Click Stops |
| Stage | Double Layer Mechanical Stage 142x132mm, Graphite Surface, Moving Range 75x50mm |
| Condenser | Abbe Condenser N.A.1.25, Center Adjustable |
| Light Source | 3W LED Brightness Adjustable |
| Adapter | C-Mount 1.0x |
| Package | 53*32.5*53.5cm/11.5Kg, 1 Set/Carton |
| 9.7" LCD Pad Digital Camera (A59.3520) | |
| LCD | 9.7" HD IPS Touch Screen With Backlight |
| Camera | 5.0M 1/2.5" CMOS Digital Camera |
| Resolution | 2048x1536 @15FPS |
| Pixel | 2.2umx2.2um |
| Data Storage | TF Card or Built-in Memory |
| Data Output | USB2.0, HDMI, TF Card |
| Photo | 5.0M .JPEG |
| Video | 1080p Video Record |
| System | Android 5.1 System |
| | 90 Languages |
| | Support Mouse + Keyboard Control |
| | Support WIFI |
| Function | Support Bluetooth |
| | S-Eye Software & Instruction Manual Installed |
| | Mark: Point, Crosshair, Coordinate, Text |
| | Length Measure: Straight Line Length, Broken Line Length, Curve Length, Parallel Line Distance, Point To Line Distance |
| | Distance Between 2 Circle Centers: Radius Fixed Circle, Two-point Fixed Circle, Three-point Fixed Circle |
| | Geometric Measure: Line Segment Length, Radius Fixed Circle, Two-point Fixed Circle, Three-point Fixed Circle, Concentric Circle |
| | Geometric Area Measure: Polygon, Square |
| Power | DC 12V/2A Power Adaptor, AC 100~240V |
| Size | Pad 238x182x12mm, Camera 40x39mm |



 0086 13911110627  sale@optoedu.com  cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China