



## Halogen Light Student 1000x WF10X Biological Compound Microscope

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1313
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Eyepiece: WF10X
- Objective: Achromatic 4x10x40x100x
- Noaeipiece: Quadruple Revolving Nosepiece
- Stage: Double Layers Mechanical Stage  
125X130mm/60X30mm
- Condenser: Abbe Condenser:N.A.1.25
- Focusing: Coarse & Fine Focusing Knob, Range 20mm
- Light Source: Halogen Lamp 6V/20W
- Highlight: **1000x Biological Compound Microscope,**  
**WF10X Biological Compound Microscope,**  
**Halogen Light wf10x microscope**



### More Images



## Product Description

OPTO-EDU A11.1313-B 1000x Binocular student microscope / Biological Compound Microscope  
Seidentopf Binocular Head, Quadruple, Halogen Lamp 6V/20w  
Seidentopf Binocular, WF10x, 25x, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 115x120mm, LED

Our Product Introduction

**OPTO-EDU**

## A11.1313 Student Biological Microscope


Binocular | Achromatic  
| LED Light



Eyepiece Diopter  
Adjustable  
Head 360°  
Rotateable



for more products please visit us on [cnoec.com](http://cnoec.com)



Precise Moving Stage | Achromatic Objective

## Abbe Condenser To Supply Bright Flat Illumination

N.A.1.25 | With Iris Diaphragm | With Filter



LED Long Life  
Light Source Adjustable



Coaxial Coarse & Fine  
Focusing Knob

## A12.1303 Series With 5 Kinds Head






A12.1303-M

A12.1303-B

A12.1303-V

A12.1303-T

A12.1303-S



A11.1313 Biological Microscope		M	S	V	B	T
Head	Monocular Head, 45° Inclined, 360° Rotatable	●				
	Compensation Free Binocular Head				●	
	Compensation Free Trinocular Head					●
	Dual Observation Head, One Tube 30° Inclined And Another Tube Vertical			●		
	Dual Observation Head, 45° Inclined		●			
Eyepiece	WF10X	●	●	●	●	●
Objective	Achromatic 4X	●	●	●	●	●
	Achromatic 10X	●	●	●	●	●
	Achromatic 40X(S)	●	●	●	●	●
	Achromatic 100X(S,Oil)	●	●	●	●	●
Nosepiece	Quadruple Revolving Nosepiece	●	●	●	●	●
Stage	Double Layers Mechanical Stage 125X130mm/60X30mm	●	●	●	●	●
Condenser	Abbe Condenser:N.A.1.25, Φ2-Φ30mm Iris Diaphragm And Φ32 Filter	●	●	●	●	●
Focusing	Coarse & Fine Focusing Knob, Range 20mm	●	●	●	●	●
Light Source	Halogen Lamp 6V/20W	●	●	●	●	●
Optional Accessories						
Light Source	LED Lamp1W					




**Opto-Edu (Beijing) Co., Ltd.**

---

 0086 13911110627

 [sale@optoedu.com](mailto:sale@optoedu.com)

 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China