



## 9" LCD Digital Zoom Stereo Optical Microscope

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A23.3645N-R1L1
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- LCD Camera: 9" LCD Screen + Digital Camera
- Head: Trinocular Head, Light Split Fixed 50:50
- Eyepiece: WF10x/20mm
- Zoom Ratio: 1:6.5
- Zoom Lens: 0.7~4.5x
- Magnification: 7x-45x
- Working Distance: 100mm, With Auxiliary Objectives 30~165mm
- Field Of View: 1~95.2mm
- Stand: Pole Stand, Size 200\*255\*22mm
- Highlight: Zoom Stereo Optical Microscope, LCD Stereo Optical Microscope, 9" Stereo Optical Microscope



### More Images



## Product Description

### OPTO-EDU A23.3645N-R1L1 7x-45x Trinocular Zoom Stereoscopic Mobile Microscope

A23.3645N-R1 Series DIY Phone PCB Repair Microscope		A23.3645N-R1	A23.3645N-R1L1	A23.3645N-R1H1
Head	Trinocular Head, With Eyepiece Adapter	•		
	LCD Screen 9" With Digital Camera, AV Output		•	
	HDMI+USB Digital Camera 38.0M Pixels			•
	45° Inclined, 360°Rotatable, Both Tube ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm			
Eyepiece	WF10x/20mm, Both With Eye-cup			
Zoom Lens	0.7~4.5x			
Zoom Ratio	1:6.5			
Magnification	7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Lens			
W.D.	Working Distance 100mm, With Standard 1x Lens			
F.O.V	Field Of View 1~95.2mm			
Stand	Pole Stand, Size 200*255*22mm, Pole Height 240mm, Dia.32mm, With Heat Resist Soft Rubber Pad, Working Area For Phone PCB Repair			
Light Source	Upper 56 LEDs Right Light Brightness Adjustable			
Optional Accessories			Cata.No.	
Eyepiece	WF15x/15mm, High-eyepoint		A51.3621-1515	
	WF20x/10mm, High-eyepoint		A51.3621-2010	
	WF25x/9mm, High-eyepoint		A51.3621-2509	
	WF30x/8mm, High-eyepoint		A51.3621-3008	
Auxiliary Lens	0.5X,W.D. 165mm		A52.3621-0.5x	
	1.5X,W.D. 45mm		A52.3621-1.5x	
	2.0X,W.D. 30mm		A52.3621-2.0x	

Objective	Standard 1.0x		Auxiliary 0.5X		Auxiliary 1.5X		Auxiliary 2X	
W,D.	100mm		165mm		45mm		30mm	
Eyepiece	Mag	View Field	Mag	View Field	Mag	View Field	Mag	View Field
10X/20mm	7.0X	28.6mm	3.5X	57.1mm	10.5X	19mm	14X	14.3mm
	45.0X	4.4mm	22.5X	8.9mm	67.5X	3mm	90X	2.2mm
15X/15mm	10.5X	21.1mm	5.25X	42.8mm	15.75X	14.3mm	21X	10.7mm
	67.5X	3.3mm	33.75X	6.7mm	101.25X	2.2mm	135X	1.7mm
20X/10mm	14.0X	14.3mm	7.0X	28.6mm	21.0X	9.5mm	28X	7.1mm
	90.0X	2.2mm	45.0X	4.4mm	135.0X	1.5mm	180X	1.1mm



Upgrade With HDMI  
Or LCD Digital  
Camera



Featuring Simul-Focal Head,  
View Image In Eyepiece &  
Camera At Same Time

## All Kinds Models Optional



A23.3645N-R1



A23.3645N-R1H1



A23.3645N-R1L1





56 LEDs Ring Light  
Brightness Adjustable


BACK HOME 



Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China