



Light 40x Eyepiece Lens Monocular Compound Optical Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1321
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Head: Monocular Head, 360° Rotatable
- Eyepiece: WF10x
- Objective: Achromatic 4x, 10x, 40x
- Nosepiece: Backward Triple Revolving Nosepiece
- Stage: Clip Stage With Disc Diaphragm
- Focusing: Coaxial Coarse & Fine Focusing Knob
- Illumination: LED Light
- Power Supply: External Transformer, 3xAA Battery Available
- Highlight: **40x compound optical microscope, eyepiece lens compound optical microscope, WF10x monocular light microscope**



More Images



Product Description

Student LED Electron Monocular Biological Microscope

Monocular, WF10x, Triple, Achromatic 4x10x40, Clip Stage 95x95mm, Single Lens Condenser + Disc Diaphragm, Coarse & Fine Focusing, LED Light Brightness Adjustable, External AD Adapter (or 3xAA Battery)
Classical & Popular Student Microscope Model In Amazon Market
Monocular Head 45° Inclined, 360° Rotatable
Eyepiece WF10x/18mm, Achromatic Objective 4x,10x,40x Spring
Plan Stage With Clips 90x90mm
LED Light For Bottom Transmit Illumination

THE KEY TO MICRO WORLD

OPTO-EDU



A11.1321 Student Biological Microscope



WF10X Eyepiece With Pointer




Achromatic Objectives 4x, 10x, 40x


Product Details



Backward triple revolving nosepiece
Achromatic 4x,10x,40x



Single layer working stage with slide clip



Disc diaphragm, adjust the amount of incoming
light to obtain higher quality images



Built-in Bottom LED light, brightness adjustable

A11.1321 Student Biological Microscope

Head	Monocular Head, 360° Rotatable
Eyepiece	WF10x
Objective	Achromatic 4x,10x,40x
Nosepiece	Backward Triple Revolving Nosepiece
Stage	Clip Stage With Disc Diaphragm
Focusing	Coaxial Coarse & Fine Focusing Knob
Illumination	LED Light
Power Supply	External Transformer, 3xAA Battery Available

A13.1321 Display






Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China