



Dual Head WF12.5x Operating Microscope For Orthopedics

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A41.1940
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Main Usage: Orthopedics, Hand, Chest, Burn And Plastic Surgery, Andrology, Genitourinary, Animal Anatomy
- 2 Heads 2 People: Main Head Binocular, 45° Inclined, Interpupillary Distance 53~76mm
- Diopter: Both Left & Right Eyetube Diopter Adjustable $\pm 6D$
- Eyepiece: WF12.5x
- Zoom Body: Motorized Continues Zoom 1:6, Parallel Light System
- Total Magnification: 2x~30x
- Highlight: **dual head operating microscope , operating microscope for orthopedics, WF12.5x operating microscope**



More Images



Product Description

Operating Microscope, Dual Head, For Orthopedics, Hand, Chest, Burn And Plastic Surgery, Andrology, Genitourinary, Animal Anatomy

A41.1940 Series Operating Microscopes, is advanced high quality series. It use fibre-optic cable for coaxial illumination, with 3 steps zoom body, has total magnification 2x~30x, featured with Motorized Focusing, it is widely used in hospital for Orthopedics, Hand, Chest, Burn And Plastic Surgery, Andrology, Genitourinary, Animal Anatomy and etc.

- 2 People 2 Binocular Head Same Magnification Same View Field
- Manual 4 Steps Zoom Body 0.6x,1.0x, 1.0x, 1.6x, Parallel Light System
- 8 Directions Motorized Foot Switch Control Focusing Up/Down, X/Y Moving And Zoom In/Out
- Four Feet With Locking Wheels, Max Extended Length 1250mm
- Coaxial Illumination, Light Spot Dia.50mm, Max Brightness 100000 Lux



OPTO-EDU (BEIJING) CO., LTD.



F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A41.1940

**Operating Microscope, Dual Head, 2x~30x
For Orthopedics, Hand, Chest, Burn And Plastic Surgery,
Andrology, Genitourinary, Animal Anatomy**





Multi Coatings Lens
Multi-layer coating, apochromatic design, high-quality image resolution, contrast and depth of field



Wide Field Eyepieces
With high eyepoint diopter adjustment function, meet the requirements of doctors wearing glasses



Continuous Magnification
Electric several manual stepless zoom, the magnification can be continuously adjusted from 2-30x

A41.1940 Surgical Operation Microscope Specification	
Main Usage	Orthopedics, Hand, Chest, Burn And Plastic Surgery, Andrology, Genitourinary, Animal Anatomy
2 Heads 2 People	Main Head Binocular, 45° Inclined, Interpupillary Distance 53~76mm Assitant Head Binocular, 45° Inclined, Interpupillary Distance 53~76mm, 180° Oppsite To Main Head, Same Magnification Same View Field As Main Head
Diopter	Both Left & Right Eyetube Diopter Adjustable $\pm 6D$
Eyepiece	WF12.5x
Zoom Body	Motorized Continues Zoom 1:6, Parallel Light System
Objective	F=200mm
Total Magnification	2x~30x
View Field	Dia.55mm~6mm
Focusing	Motorized Focusing Controlled By Foot Switch, Focusing Range +/-25mm, Focusing Speed 33s/Full Range, Auto Return To Zero
Foot Switch	8 Directions X/Y/Z + Zoom In/Out Motorized Foot Switch --Control Focusing Up/Down, Moving Range 50mm --Control X/Y Moving, Moving Range 60mm --Control Magnification Zoom In/Out --Water Proof Design

Stand	Luxury Stand With Control Electrical Parts Built-in --Max Extended Length 1250mm --Max Rotatable Angle Of Arm 0~300° --Distance Between Objective And Gound Surface 800mm~1250mm
Illumination	Coaxial Illumination, Light Spot Dia.50mm, Max Brightness 100000 Lux
Weight	164Kgs



A52.1941 Objective



A52.1941-1830



A55.1941



A55.1942

Optional Accessories		Cata. No.
Objective	F=175mm	A52.1941-175
	F=250mm	A52.1941-250
	F=300mm	A52.1941-300
	F=350mm	A52.1941-350
	Zoom F=180~300mm	A52.1941-1830
Adapter	Beamsplitter	A55.1941
	C-Mount 0.3x For 1/3" Digital Camera	A55.1942



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China