



Stereo OPTO EDU 20x-100x Binocular Biological Microscope

Our Product Introduction


Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1526
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month




Product Specification

- Magnification: 20x~100x
- Focusing: Coarse Focusing Knob, Focusing Range 25mm
- Power: 3xAAA Battery
- Packing: Polyfoam Protection + Strong Master Carton
- Highlight: **100x Binocular Biological Microscope, 20x Binocular Biological Microscope, OPTO EDU Binocular Biological Microscope**



More Images



for more products please visit us on cnoec.com

Product Description

Biological & Stereo Microscope, Monocular, WF10x/18mm, Zoom Lens 2-10x, Round Stage Dia.82mm, Disc Diaphragm, Up/Bottom LED, Cordless 3xAAA Batteries

Our Product Introduction


THE KEY TO MICRO WORLD

OPTO-EDU





A11.1526

Student Biological & Stereo Microscope



Biological & Stereo
Observation All-in-One

Up & Bottom LED,
Brightness Adjustable



Cordless Power
Supply By 3xAAA
Batteries



A11.1526 Student Biological & Stereo Microscope	
Magnification	20x~100x
Focusing	Coarse Focusing Knob, Focusing Range 25mm
Power	3xAAA Battery
Packing	Polyfoam Protection + Strong Master Carton




for more products please visit us on cnoec.com

Product Details




Monocular Head, With WF10x/18mm Eyepiece



Zoom 2-10x Objective, Upper 10 LEDs Light



Bottom Illumination LED, Brightness Adjustable



Round Clip Stage Dia.82mm With Iris Diaphragm



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China