



## Opto Edu 40x-640x Monocular Compound Microscope

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1527
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Head: Monocular, 45°
- Eyepiece: WF10x-16x Zoom Eyepiece
- Nosepiece: Triple
- Objective: 4x, 10x, 40xs
- Magnification: 40x~640x
- Working Stage: Clip Stage, 81x81mm, With 6 Holes Iris Diaphragm
- Focusing: Coarse Focusing Knob, Focusing Range 25mm
- Light Source: LED Illumination, Brightness Adjustable
- Power: 3xAAA Battery
- Packing: Polyfoam Protection + Strong Master Carton
- Highlight: **opto edu monocular compound microscope ,  
40x monocular compound microscope ,  
640x monocular compound microscope**



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)


## Product Description

Monocular, 45°, WF10x-16x zoom eyepiece, Triple, Objective 4x10x40xs, stage 81\*81mm with clip, 6 holes disc disphragm, LED lamp, brightness adjustable, plastic base with battery pack.

A11.1527 Student Biological Microscope	
Head	Monocular, 45°
Eyepiece	WF10x-16x zoom eyepiece
Nosepiece	Triple
Objective	4x, 10x, 40xs
Magnification	40x~640x
Working Stage	Clip Stage, 81x81mm, With 6 Holes Iris Diaphragm
Focusing	Coarse Focusing Knob, Focusing Range 25mm
Light Source	LED Illumination, Brightness Adjustable
Power	3xAAA Battery
Packing	Polyfoam Protection + Strong Master Carton

THE KEY TO MICRO WORLD

OPTO-EDU




## A11.1527 Student Biological Microscope



Triple Nosepiece,  
Objective 4x10x40xs

Cordless LED Lamp,  
Brightness Adjustable



Plastic Base With  
Battery Pack



Product Multi-angle



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China