



OPTO EDU A11.1009-E Sliding Binocular Compound Microscope For Student

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1009
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Head: Sliding Binocular Head, 45° Inclined, 360° Rotatable
- Eyepiece: Wide Field WF10x/18mm
- Nosepiece: Quadruple Nosepiece
- Achromatic Objective: 4x/0.10, 10x/0.25, 40x/0.65(S), 100x/1.25(S,Oil)
- Magnification: 40x~1000x
- Stage: Double Layer Mechanical Stage, 110*120mm, Moving Range 70*20mm
- Focusing: Coaxial Coarse & Fine Focusing
- Condenser: Abbe Condenser N.A.1.2, With Iris Diaphragm, With Blue Filter
- Light Source: LED Lamp Illumination 1W, Brightness Adjustable
- Spare Parts: Dust Cover, Spare Bulb & Fuse
- Highlight: **OPTO EDU binocular compound microscope , Sliding binocular light microscope ,**



Product Description

High Quality Compact Design Student Biological Microscope
Monocular, Binocular Head 45° Inclined, 360° Rotatable, WF10x/18mm Eyepiece
3/4 Holes Nosepiece With Achromatic Objective 4x10x40x100x
Mechanical Stage, Abbe Condenser, Coaxial Coarse & Fine Focusing
Coaxial Coarse & Fine Focusing Knob, Light Source Brightness Adjustable



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A11.1009

Student Biological Microscope



Our Product Introduction

for more products please visit us on cnoec.com



A11.1009-D Monocular Head



A11.1009-E Binocular Head




Achromatic Objectives



Double Layer Mechanical Stage



Brightness Adjustable



Abbe Condenser N.A.1.2, With Iris Diaphragm



A11.1009-D



A11.1009-E

Items	A11.1009 Student Biological Microscope	M	B
Head	Monocular Head, 45° Inclined, 360° Rotatable	●	
	Sliding Binocular Head, 45° Inclined, 360° Rotatable		●
Eyepiece	Wide Field WF10x/18mm	●	●
	WF5x	○	○
	WF16x	○	○
	WF10x Eyepiece With Pointer	○	○
Nosepiece	Triple Nosepiece	●	
	Quadruple Nosepiece		●
Achromatic Objective	4x/0.10, 10x/0.25, 40x/0.65(S)	●	
	4x/0.10, 10x/0.25, 40x/0.65(S), 100x/1.25(S,Oil)		●
	20x 60x	○	○
Magnification	40x~400x	●	
	40x~1000x		●
Stage	Double Layer Mechanical Stage, 110*120mm, Moving Range 70*20mm	●	●
Focusing	Coaxial Coarse & Fine Focusing	●	●
Condenser	Abbe Condenser N.A.1.2, With Iris Diaphragm, With Blue Filter, Rack & Pinion Up/Down Adjustment	●	●
Light Source	Incandescent Lamp Illumination 15W, Brightness Adjustable	●	
	LED Lamp Illumination 1W, Brightness Adjustable		●
	Mirror	○	○
Spare Parts	Dust Cover, Spare Bulb & Fuse	●	●

Note: "●" In Table Is Standard Outfits, "○" Is Optional Accessories.



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China