



## Monocular Triple Wf10x Microscope With Rechargeable Battery

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1133
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



### Product Specification

- Head: Monocular Head, 30° Inclined, 360° Rotatable
- Eyepiece: WF10X With Pointer
- Objective: Achromatic 4X 10X 40X
- Nosepiece: Triple Revolving Nosepiece
- Stage: Plain Stage With Clips, Size 100x90mm
- Condenser: Single Lens Condenser With Disc Diaphragm
- Focusing: Coaxial Coarse And Fine Focusing Knobs
- Light Source: LED Lampe With Rechargeable Battery
- Highlight: **triple wf10x microscope, monocular wf10x microscope, wf10x microscope with rechargeable battery**

OPTO-EDU



### More Images

OPTO-EDU



OPTO-EDU



## Product Description

### Monocular Triple Wf10x Microscope With Rechargeable Battery

Items	A11.1133 Biological Microscope Specification	A	B	C	D
Head	Monocular Head, 30° Inclined, 360° Rotatable	●	●		
	Sliding Binocular Head, 45° Inclined, 360° Rotatable			●	
	Seidentopf Binocular Head, 30° Inclined, 360° Rotatable				●
Eyepiece	WF10X With Pointer	●	●	●	●
	Achromatic 4X	●	●	●	●
Objective	Achromatic 10X	●	●	●	●
	Achromatic 40X, Spring	●	●	●	●
	Achromatic 100X, Spring/Oil	●	●	●	●
	Triple Revolving Nosepiece	●			
Stage	Quadruple Revolving Nosepiece	●	●	●	●
	Plain Stage With Clips, Size 100x90mm	●			
	Attachable Mechanical Stage, Size 100x90mm		●	●	
Condenser	Double Layer Mechanical Stage, Size 100x90mm				●
	Single Lens Condenser With Disc Diaphragm	●			
Light Source	N.A.1.25 Condensor With Iris Diaphragm		●	●	●
	Coaxial Coarse And Fine Focusing Knobs	●	●	●	●
Light Source	LED Lampe With Rechargeable Battery	●	●	●	●
	Power Supply Adapter	●	●	●	●



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China  
Tel:+8610 88696020 Fax:+8610 88696085

## A11.1133 Biological Microscope



Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Display



A11.1133-A



A11.1133-B



A11.1133-C



A11.1133-D

## Product Specification



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China