



## 40x Student Led Electron Monocular Biological Microscope

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1321
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Head: Monocular Head, 360° Rotatable
- Eyepiece: WF10X
- Objective: Achromatic 4x,10x,40x
- Nosepiece: Backward Triple Revolving Nosepiece
- Stage: Clip Stage With Disc Diaphragm
- Focusing: Coaxial Coarse & Fine Focusing Knob
- Illumination: LED Light
- Power Supply: External Transformer, 3xAA Battery Available
- Highlight: **electron monocular biological microscope, led monocular biological microscope, 40x monocular biological microscope**



### More Images




## Product Description

|  |  |
|--|--|
| A11.1321 Student Biological Microscope |  |
| Head                                   | Monocular Head, 360° Rotatable               |
| Eyeiece                                | WF10x  |
| Objective                              | Achromatic 4x,10x,40x                        |
| Nosepiece                              | Backward Triple Revolving Nosepiece          |
| Stage                                  | Clip Stage With Disc Diaphragm               |
| Focusing                               | Coaxial Coarse & Fine Focusing Knob          |
| Illumination                           | LED Light                                    |
| Power Supply                           | External Transformer, 3xAA Battery Available |


## Product Details



Backward triple revolving nosepiece  
Achromatic 4x, 10x, 40x



Disc diaphragm, adjust the amount of incoming light to obtain higher quality images



Single layer working stage with slide clip



Built-in Bottom LED light, brightness adjustable



WF10X Eyepiece With Pointer



Achromatic Objectives 4x, 10x, 40x

## Product Display



## Product Specification



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627    ✉ sale@optoedu.com    🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China