



Wifi 2.0m Android Digital Biological Microscope With 6-13" Pad Holder A31.0921-2M

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: CHINA
- Brand Name: OPTO-EDU
- Certification: CE, ISO
- Model Number: A31.0921-2M
- Minimum Order Quantity: 1
- Price: FOB \$1~1000, Depend on Order Quantity
- Delivery Time: 15-25days
- Payment Terms: T/T, L/C, D/A, D/P, Western Union, MoneyGram
- Supply Ability: 100pcs per month



Product Specification

- Magnification: 40x~400x
- Head: Monocular Head, 30° Inclined Built-in WIFI Camera, 2MP. With Pad Holder
- Eyepiece: PL10X/18mm
- Nosepiece: Triple Nosepiece
- Objective: Achromatic 4x/0.10, 10x/0.25, 40xs/0.65
- Stage: Metal Plain Stage 110*120mm, With Slide Clip
- Focusing: Separate Coarse & Fine Focusing. Coarse Focusing Range 8.7mm, Fine Focusing Range 1.2mm
- Condenser: NA0.65 Condenser And Iris Diaphragm
- Light Source: LED Illumination (3 Bulbs)
- Highlight: **wifi digital biological microscope, 2.0m digital biological microscope, android digital biological microscope**



Product Description

A31.0921	
Magnification	40x~400x
Head	Monocular Head, 30° Inclined Built-in WIFI camera, 2MP. With Pad holder
Eyepiece	PL10X/18mm
Nosepiece	Triple Nosepiece
Objective	Achromatic 4x/0.10, 10x/0.25, 40xs/0.65
Stage	Metal Plain Stage 110*120mm, With Slide Clip
Focusing	Separate Coarse & Fine Focusing. Coarse Focusing Range 8.7mm, Fine Focusing Range 1.2mm
Condenser	NA0.65 condenser and iris diaphragm
Light Source	LED Illumination (3 bulbs)



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A31.0921-2M WIFI Digital Microscope



Product Details



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China