



## A12.1010 OPTO EDU 40x Biological Microscope 20mm Eyepiece

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE
- Model Number: A12.1010
- Minimum Order Quantity: 1
- Price: \$1-\$1000
- Packaging Details: Carton
- Delivery Time: 1month
- Payment Terms: T/T, West Union and Paypal
- Supply Ability: 10000



### Product Specification

- Eyepiece: Wide Field WF10x/20mm
- Nosepiece: Backward Quadruple Nosepiece
- Condenser: N.A.1.25 Abbe Condense, Iris Diaphragm
- Working Stage: Double Layer Mechanical Stage, 216x150mm
- Highlight: **laboratory opto edu microscope, 20mm eyepiece biological microscope, biological microscope 40x**



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

Binocular & Trinocular, WF10x/20, Quadruple, Infinity Plan 4x10x40x100x, Mechanical Stage, Abbe 1.25, Coaxial Coarse&Fine, 12V20W

Items	A12.1010 Laboratory Biological Microscope	B	T
Head	Seidentopf Binocular Head, 30° Inclined, 360°Rotatable, Interpupillary Distance 48-75mm	•	
	Seidentopf Trinocular Head, 30° Inclined, 360°Rotatable, Interpupillary Distance 48-75mm		•
Eyepiece	Wide Field WF10x/20mm, Eyepiece Tube Dia.30mm, 2 pcs	•	•
Nosepiece	Backward Quadruple Nosepiece	•	•
Optical System	Infinity Optical System	•	•
Objective	Infinity Plan Achromatic 4x/0.10, 10x/0.25, 40x/0.65(S), 100x/1.25(S, Oil)	•	•
Magification	40x-1000x	•	•
Stage	Double Layer Mechanical Stage, 216x150mm, Moving range 54x78mm, Double Slide Holder, Low Position X/Y Coaxial Control Knob, With Travel Guide To Avoid Scratch	•	•
Focusing	Low Position Coaxial Coarse Fine Focusing.	•	•
	Fine Focusing Scale 0.002mm, Moving Range 25mm, With Safety Stop	•	•
Condenser	N.A.1.25 Abbe Condense, Iris Diaphragm, With Filter Holder, Up/Down Range 25mm, Center Adjustable, With Filter Holder	•	•
Filter	Green	•	•
Light Source	Halogen Lamp 12V/20W, Brightness Adjustable	•	•
Note: "•" In Table Is Standard Outfits			
Items A12.1010 Laboratory Biological Microscope Optional Accessories			
Eyepiece	Wide Field WF15x/16mm, Eyepiece Tube Dia.30mm	A51.1002-15B	
	Wide Field WF20x/12mm, Eyepiece Tube Dia.30mm	A51.1002-20B	
Objective	Infinity Plan Achromatic 20x	A52.1008-20	
	Infinity Plan Achromatic 60x	A52.1008-60	
Nosepiece	Backward Quintuple Nosepiece	A54.1001-5	
Illumination	Kohler Illumination	A56.1005-5	
Adapter	Photo Adapter	A55.1008	
	C-Mount 0.4x	A55.1009-04	
	C-Mount 0.5x	A55.1009-05	
Dark Field	Dark Field Condense, Dry	A5D.1012-D	
	Dark Field Condense, Immersion	A5D.1012-I	
Polarizing	Simple Polarizing Kit	A5P.1005	
Phase Contrast	Phase Contrast Set	A5C.1006	
Fluorescent	EPI-Fluorescent Kit	A5F.1005	



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China  
Tel: +8610 88696020 Fax: +8610 88696085

# A12.1010


## Laboratory Biological Microscope



## A12.1010 Laboratory Biological Microscope



A12.1010-A  
Binocular, Quadruple  
Laboratory Biological Microscope




A12.1010-B  
Trinocular, Quadruple  
Laboratory Biological Microscope



Seidentopf Trinocular Head, 30° Inclined, 360° Rotatable, Interpupillary 48~75mm, Wide Field WF10x/20mm, Dia.30mm



Backward Quadruple Nosepiece, With Infinity Plan Achromatic 4x/0.10, 10x/0.25, 40x/0.65(S), 100x/1.25(S,Oil)



Double Layer Mechanical Stage, 216x150mm, Moving Range 54x78mm, Double Slide Holder, Low Position X/Y Coaxial Control Knob, With Travel Guide To Avoid Scratch



Halogen Lamp 12V/20W, Brightness Adjustable  
Light Source Cover Removable For Easy Repair





N.A.1.25 Abbe Condense, Iris Diaphragm, With Filter Holder, Up/Down Range 25mm, Center Adjustable, With Filter Holder



You Can See The Change Of Objective Magnification Even From The Back



**Focusing**  
Low Position Coaxial Coarse & Fine  
Focusing




**Knob**  
Switch & Light Source Adjustment knob



**Opto-Edu (Beijing) Co., Ltd.**

 0086 13911110627

 [sale@optoedu.com](mailto:sale@optoedu.com)

 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China