



## 2.0m 1.3m 7" Digital Lcd Microscope Stereo Usb Portable Dual Lens

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)


### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A36.5101
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- LCD Screen: 7" IPS Screen, With Touch Button To Take Photo & Video
- Handheld Camera: 1.3M Digital Camra, Zoom Lens 1.5x~25x
- Main Body Camera: 2.0M Digital Camra, Zoom Lens 5x~35x
- TF Card Output: Support Max 32G, To Take Photo & Video
- USB Output: To Win/Mac, With Professional StrangeView Measuring Software
- Stand: Pole Stand, With Focusing Knob, Size 180x170mm
- Size: 180\*170\*322mm
- Power Supply: DC 5V2A Adapter, Can Be Powered By USB Power Bank For Outdoor Use
- Highlight: **2.0m digital lcd microscope, 7" stereo digital lcd microscope, 7 inch lcd digital usb microscope**



### More Images



## Product Description

### 7" LCD Stereo + USB Portable Dual Lens Digital Microscope, 2.0M+1.3M

Brand New Design Dual Lens LCD + Portable Digital Microscope For Different Purposes

7.0" LCD, Main Body 2.0M Camera, 5x-35x Zoom For Real Time Clear Watching

Portable USB2.0 1.3M Digital Microscope, 1.5x-25x Zoom For Large Size Objective Watch on LCD

USB Output To Win/Mac, With Professional StrangeView Measuring Software

Support Max 32G TF Card To Take Photo & Video, Power By 5V2A Adapter



Metal fuselage



Plated film lens



Taking pictures



Connect the computer



LED light source



Touch control

**HD**

16 million pixels



Adjustable zoom



Super clear quality

## Product Details



## Metal Zoom Lens 5x~35x

Aviation Aluminum, CNC Carved Molding  
High-definition Coating Zoom Lens, The  
Focus Is More Precise, The Image Is Clearer



## Handheld Viewing

Is More Convenient  
Large Size PCB Can Also Be Tested  
Quickly



## View on Computer

Good Helper For Electronic Repair  
Dual-lens Design, More Convenient For  
Detection And Storage



## LED Illumination


Stepless Dimming Of The Main  
Lens Hold The Lens And Touch The  
Button To Dim






## Metal body



Aviation Aluminum, Sturdy And Durable



## Practical Application



Clock & Electronic Precision Repair





Lapis Lazuli



Emerald Rough Stone

Gem Detection




Detection Of Various Currencies & Commemorative Coins & Banknotes & Stamps




Textile Testing

A36.5101 7" LCD Stereo + USB Portable Dual Lens Digital Microscope, 2.0M+1.3M	
LCD Screen	7" IPS Screen, With Touch Button To Take Photo & Video, One Key Switch Between Main Body Camera & Portable USB Camera
Handheld Camera	1.3M Digital Camra, Zoom Lens 1.5x~25x, Output to 7" LCD Screen 8 LEDs Light, Brightness Adjustable
Main Body Camera	2.0M Digital Camra, Zoom Lens 5x~35x, TF Card, USB Output Coaxial LED Ring Light, Brightness Adjustable
TF Card Output	Support Max 32G, To Take Photo & Video
USB Output	To Win/Mac, With Professional StrangeView Measuring Software
Stand	Pole Stand, With Focusing Knob, Size 180x170mm
Size	180*170*322mm
Power Supply	DC 5V2A Adapter, Can Be Powered By USB Power Bank For Outdoor Use
Package	Color Box + Foam, 1 pc/1 Box, With Software Disc, DC Adapter, Instruction Manual

PC Measuring Software



## LCD Screen Touch Menu



## Product Packing






**Opto-Edu (Beijing) Co., Ltd.**

---

 0086 13911110627

 [sale@optoedu.com](mailto:sale@optoedu.com)

 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China